Cavitron® Plus Ultrasonic Scaler

with Tap-On™ Technology

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Please read carefully and completely before operating unit.

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

Quick Start Installation Instructions


Cavitron®


- 1** Set unit to **OFF** position.
Symbol: ○ means off.


- 2** Connect **POWER** supply.


- 3** Connect **WATER** supply by pushing **BLUE** water hose into water port.


- 4** Connect **AIR** supply (if equipped) by pushing **BLACK** air hose into air port.

- 5** Install **BATTERIES** in wireless foot pedal unit.


- 6** Connect **HANDPIECE** to cable.
TIP: Handpiece is detachable and autoclavable. Make sure the electrical connections are dry before connecting.


- 7** Set unit to **ON** position.
Symbol: | means on.
The indicator light will be green when unit is on.


- 8** Hold upright and depress foot pedal until handpiece **FILLS WITH WATER**.
TIP: adjust lavage control to higher flow rate to increase water flow.
Tap-On™ Technology does not function unless an insert is in the handpiece.

- 9** Place **ULTRASONIC INSERT** into handpiece.

- 10** Adjust **POWER LEVEL**.

- 11** Adjust **LAUGE CONTROL** on handpiece to desired setting.


- 12** Consult **Directions For Use** for daily maintenance, usage techniques and more.


Quick Start User Guide

Cavitron®


Tap-On™ Technology

To **activate** Tap-On™ Technology:
Tap wireless foot pedal one time quickly.

Relax foot as you scale or air polish.

To **deactivate** Tap-On™ Technology:
Tap wireless foot pedal one time.

User controls

Turbo: When pressed, increases power delivered to the system by up to 25%; turbo power remains on until the button is pressed again.

Boost: Temporary, hands-free activation by pressing the foot pedal all the way to the floor; enables quick removal of tenacious calculus.

Tip: Quick release of pedal from Boost will not deactivate Tap-On™ Technology.

Purge: Remove insert from Handpiece and press Purge button; water will purge through system for two minutes or until you press the foot pedal or press the Purge button again.

Finer Lavage Water Control:

Conveniently adjust lavage water control directly on the handpiece;

- 1 = lowest water flow
- 6 = highest water flow
- > 6 = flush or purge

Quick Reference Guide Diagnostic Display

ON/OFF

Illuminates when the Main Power On/Off Power switch is in the “ON” (I) position.

TURBO

Offers the ability to increase power to the system by up to 25% with the push of a button. Purple arrows illuminate when in use.

BOOST

Illuminates when the Boost Mode is activated by the Tap-On™ Wireless Foot Pedal. To activate, fully depress Tap-On™ Foot Pedal to the second position (all the way to the floor). To deactivate, release Tap-On™ Foot Pedal to the first position.

PURGE BUTTON

Illuminates when Purge function is activated. To activate Purge, remove insert from handpiece, turn the Handpiece Lavage Control to maximum water flow, press the Purge Button on the Diagnostic Display. Water will purge through system lines for two minutes. To deactivate during two minute cycle, press Purge button again or press Tap-On™ Foot Pedal.

SERVICE

Illuminates when the system is not functioning properly. This display has three distinct modes:


- Slow blink (1 blink per second) means the system is not operating within factory specifications.
- Fast blink (3 blinks per second) indicates an improper set-up.
- Steady light indicates the system is overheating.

LOW BATTERY

Illuminates when the Tap-On™ Foot Pedal battery power is approaching end of life. Replace batteries as instructed in section 7.9.

Power Control

Power Level Control

Turn knob to select ultrasonic power level for operation. Turning the knob clockwise increases the distance the insert tip moves (the stroke) without changing frequency; turning knob counter-clockwise decreases the distance the insert tip moves (the stroke) without changing the frequency.

RINSE

Rinse mode is used during an ultrasonic scaling procedure when lavage is required to flush the procedural area. To activate, turn Power Level Control Knob fully counter-clockwise until a “click” is heard.

BLUE ZONE

Provides an extended low-power range for improved patient comfort when subgingivally scaling.

TABLE OF CONTENTS

QUICK START INSTALLATION INSTRUCTIONS 2

QUICK START USER GUIDE 3

QUICK REFERENCE GUIDE: DIAGNOSTIC DISPLAY 4

INTRODUCTION 6

PRODUCT OVERVIEW 6

TECHNICAL SUPPORT 6

SUPPLIES & REPLACEMENT PARTS 6

INDICATIONS FOR USE 7

CONTRAINDICATIONS 7

WARNINGS 7

PRECAUTIONS

 4.1 System Precautions 7-8

 4.2 Procedural Precautions 8

ADVERSE REACTIONS 8

INFECTION CONTROL

 6.1 General Infection Control 8

 6.2 Water Supply Recommendations 8

INSTALLATION INSTRUCTIONS

 7.1 Water Line Requirements 8-9

 7.2 Electrical Requirements 9

 7.3 Unpacking the System 9

 7.4 System Installation 9

 7.5 Power Cord Connection 9

 7.6 Water Supply Line Connection 9

 7.7 Tap-On™ Foot Pedal Battery Installation/ Replacement 10

 7.8 Tap-On™ Foot Pedal Synchronization 10

CAVITRON® Plus Scaler DESCRIPTION

 8.1 System Controls 11

 8.2 Diagnostic Display Indicators and Controls 12

 8.3 Handpiece/Cable 13

 8.4 Cavitron® 30K™ Ultrasonic Inserts 13

 8.5 Tap-On™ Wireless Foot Pedal Operation 14

 8.6 Accessories and User Replaceable Parts 14

 8.6.1 Accessories 14

 8.6.2 User Replaceable Part Kits 14

SYSTEM SETUP, OPERATION AND TECHNIQUES FOR USE

 9.1 Handpiece Setup 14-15

 9.2 Turbo Mode 15

 9.3 Boost Mode 15

 9.4 Patient Positioning 15

 9.5 Performing Ultrasonic Scaling Procedures 15

 9.6 Patient Comfort Considerations 15

SYSTEM CARE

 10.1 Daily Maintenance 16

 Start-up procedures at the beginning of the day 16

 Between patients 16

 Shut-down procedures at the end of the day 16

 10.2 Weekly Maintenance 16

 10.3 Monthly Maintenance 17

 Water Line Filter Maintenance 17

TROUBLESHOOTING

 11.1 Troubleshooting Guide 17-18

 11.2 Technical Support and Repairs 18

WARRANTY PERIOD 18

SPECIFICATIONS 18-19

SYMBOL IDENTIFICATION 19

CLASSIFICATIONS 19

DISPOSAL OF UNIT 19

ELECTRO-MAGNETIC COMPATIBILITY PRECAUTIONS 20-22

QUICK REFERENCE GUIDE: TROUBLESHOOTING 23

INTRODUCTION

Congratulations!

Your decision to add the Cavitron® Plus Ultrasonic Scaler with Tap-On™ Technology to your practice represents a wise investment in good dentistry.

For over four decades, dental professionals have preferred the clinical benefits and labor-saving advantages inherent in Cavitron ultrasonic scalers. Clinical studies and independent research have confirmed the speed, efficiency and versatility of ultrasonic scaling.*

DENTSPLY Professional is an ISO 13485 registered company. All DENTSPLY Professional medical devices sold in Europe are CE marked in conformance with Council Directive 93/42/EEC.

Website: www.professional.dentsply.com

CAUTION: United States Federal Law restricts this device to sale by or on the order of, a licensed dental professional.

For dental use only.

PRODUCT OVERVIEW

The Cavitron® Plus Ultrasonic Scaler is a precision engineered and manufactured instrument. It contains controls and components for ultrasonic scaling. The system produces 30,000 strokes per second at the ultrasonic insert's working tip that when combined with the cavitation effect of the coolant lavage creates a synergistic action that is designed to "power away" even the heaviest calculus deposits while maintaining operator and patient comfort.

The Cavitron Plus Ultrasonic Scaler is equipped with the Sustained Performance System™ (SPS Technology), which offers a constant balance between scaling efficiency and patient comfort by maintaining power when the insert tip encounters tenacious deposits, allowing the clinician to effectively scale even at a decreased/lower power setting. The Cavitron Plus System has extended the SPS technology by expanding the Blue Zone range, providing finer resolution to the power settings.

Advanced features that make the Cavitron Plus a wise investment include a Tap-On™ Wireless Foot Pedal with Tap-On™ Technology, Turbo Mode, illuminated diagnostic display, rinse setting, automated purge function, Steri-Mate® detachable sterilizable handpiece, and 330° swivel handpiece cable with more precise lavage water control. These features, combined with a low power range (Blue Zone™) and hands-free Boost Mode, are designed to deliver a positive ultrasonic scaling experience for your patients while providing your practice with the quality and reliability you've come to expect from Cavitron brand ultrasonic systems.

The Cavitron Plus Ultrasonic Scaler is UL/ULC certified and approved. The Cavitron Plus Ultrasonic Scaler is classified by Underwriters Laboratories Inc. with respect to electric shock, fire, mechanical hazards in accordance with the IEC 60601 Standard. The Cavitron Plus Ultrasonic Scaler complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: 1) this device may not cause harmful interference, and 2) this device must accept any interference received, including interference that may cause undesired operation. Cavitron Plus base FCC certification/registration number: FCC ID: TF3-DPD81842; IC: 4681B-DPD81842. Cavitron Plus Tap-On™ Foot Pedal FCC certification/registration number: FCC ID: TF3-DPD81861; IC: 4681B-DPD81861. The term IC before the certification/registration number signifies that the Industry Canada technical specifications were met.

TECHNICAL SUPPORT

For technical support and repair assistance in the U.S., call the DENTSPLY Professional Cavitron CareSM Factory Certified Service at 1-800-989-8826, Monday through Friday, 8:00 A.M. to 5:00 P.M. (Eastern Time). For other areas, contact your local DENTSPLY® Professional representative.

SUPPLIES & REPLACEMENT PARTS

To order supplies or replacement parts in the U.S., contact your local DENTSPLY Professional Distributor or call 1-800-989-8826, Monday through Friday, 8:00 A.M. to 5:00 P.M. (Eastern Time). For other areas, contact your local DENTSPLY Professional Representative.

SECTION 1: Indications For Use

- All general supra and subgingival scaling applications
- Periodontal debridement for all types of periodontal diseases
- Endodontic procedures

SECTION 2: Contraindications

- Ultrasonic Systems should not be used for restorative dental procedures involving the condensation of amalgam.

SECTION 3: Warnings

- The use of High Volume Saliva Evacuation to reduce the quantity of aerosols released during treatment is highly recommended.
- Prior to beginning treatment, patients should rinse with an antimicrobial such as Chlorhexidine Gluconate 0.12%. Rinsing with an antimicrobial reduces the chance of infection and reduces the number of microorganisms released in the form of aerosols during treatment.
- It is the responsibility of the Dental Healthcare Professional to determine the appropriate uses of this product and to understand:
 - the health of each patient,
 - the dental procedures being undertaken,
 - and applicable industry and governmental agency recommendations for infection control in dental healthcare settings,
 - requirements, and regulations for safe practice of dentistry; and
 - these Directions for Use in their entirety, including Section 4 Precautions, Section 6 Infection Control, and Section 10 System Care.
- Where asepsis is required or deemed appropriate in the best professional judgment of the Dental Healthcare Professional, this product should not be used, unless used in combination with a sterile lavage kit (P/N 81340).
- Handle Cavitron insert with care. Improper handling of insert, specifically the insert tip, may result in injury and/or cross contamination.
- Failure to follow properly validated sterilization processes and approved aseptic techniques for Cavitron inserts may result in cross contamination.
- Persons fitted with cardiac pacemakers, defibrillators and other active implanted medical devices, have been cautioned that some types of electronic equipment might interfere with the operation of the device. Although no instance of interference has ever been

reported to DENTSPLY, we recommend that the handpiece and cables be kept 6 to 9 inches (15 to 23 cm) away from any device and their leads during use.

There are a variety of pacemakers and other medically implanted devices on the market. Clinicians should contact the device manufacturer or the patient's physician for specific recommendations. This unit complies with IEC 60601 Medical Device Standards.

- Insufficient water flow could result in elevated water and tip temperature. When operated at the input water temperature specified in the Water Line Requirements Section 7.1 and with sufficient water flow, the water and tip temperature should not exceed 50° C (122° F). Failure to follow recommendations for environmental operating conditions, including input water temperature, could result in injury to patients or users. If temperature is elevated, increase water flow. If temperature remains elevated, discontinue use.
- During boil-water advisories, this product should not be operated as an open water system (e.g. connected to a public water system). A Dental Healthcare Professional should disconnect the system from the central water source. The Cavitron DualSelect system can be attached to this unit and operated as a closed system until the advisory is cancelled. When the advisory is cancelled, flush all incoming waterlines from the public water system (e.g. faucets, waterlines and dental equipment) in accordance with the manufacturer's instructions for a minimum of 5 minutes.
- Per FCC Part 15.21, changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.
- This portable transmitter with its antenna complies with FCC/IC RF exposure limits for general population / uncontrolled exposure.
- This Device complies with Industry Canada License-exempt RSS standard(s). Operation is subject to the following two conditions: 1) this device may not cause interference, and 2) this device must accept any interference, including interference that may cause undesired operation of the device.

SECTION 4: Precautions

4.1 System Precautions

- Close manual shut-off valve on the dental office water supply every night before leaving the office.
- Do not place the system on or next to a radiator or other heat source. Excessive heat may damage the system's electronics. Place the system where air is free to circulate on all sides and beneath it.
- The system is portable, but must be handled with care when moving.

- Equipment flushing and dental water supply system maintenance are strongly recommended. See Section 10: System Care.
- Never operate system without fluid flowing through the handpiece.
- Always ensure that the electrical connections on the handpiece cable and the Steri-Mate® handpiece are clean and dry before assembling them for use.

4.2 Procedural Precautions

General

- As with all dental procedures, use universal precautions (i.e., wear face mask, eyewear, or face shield, gloves and protective gown).

Ultrasonics

- The Cavitron Plus unit works with Cavitron inserts as a system, and was designed and tested to deliver maximum performance for all currently available Cavitron brand ultrasonic inserts. Companies that manufacture, repair or modify inserts carry the sole responsibility for proving the efficacy and performance of their products when used as a part of this system. Users are cautioned to understand the operating limits of their insert before using in a clinical setting.
- Like bristles of a toothbrush, ultrasonic inserts “wear” with use. Inserts with just 2 mm of wear lose about 50% of their scaling efficiency. In general it is recommended that ultrasonic inserts be discarded and replaced after one year of use to maintain optimal efficiency and avoid breakage. A DENTSPLY Professional Insert Efficiency Indicator is enclosed for your use.
- If excessive wear is noted, or the insert has been bent, reshaped or otherwise damaged, discard the insert immediately.
- Ultrasonic insert tips that have been bent, damaged, or reshaped are susceptible to in-use breakage and should be discarded and replaced immediately.
- Retract the lips, cheeks and tongue to prevent contact with the insert tip whenever it is placed in the patient’s mouth.

SECTION 5: Adverse Reactions

None Known.

SECTION 6: Infection Control

6.1 General Infection Control

- For operator and patient safety, carefully practice the infection control procedures detailed in the Infection Control Information Booklet accompanying your

System. Additional booklets can be obtained by calling Customer Service at 1-800-989-8826, Monday through Friday, 8:00 A.M. to 5:00 P.M. (Eastern Time). For areas outside the U.S., contact your local DENTSPLY Professional representative.

- As with high speed handpieces and other dental devices, the combination of water and ultrasonic vibration from the Cavitron Plus Ultrasonic Scaler will create aerosols. Following the procedural guidelines in Section 9 of this manual can effectively control and minimize aerosol dispersion.

6.2 Water Supply Recommendations

- It is highly recommended that all dental water supply systems conform to applicable CDC (Centers for Disease Control and Prevention) and ADA (American Dental Association) standards, and that all recommendations be followed in terms of flushing, chemical flushing, and general infection control procedures. See Sections 7.1 and 10.
- As a medical device, this product must be installed in accordance with applicable local, regional, and national regulations, including guidelines for water quality (e.g. drinking water). As an open water system, such regulation may require this device to be connected to a centralized water control device. The Cavitron® DualSelect™ Dispensing System may be installed to allow this unit to operate as a closed water system.

SECTION 7: Installation Instructions

Anyone installing a Cavitron Plus System should observe the following requirements and recommendations.

7.1 Water Line Requirements


- A water supply line with user-replaceable filter is supplied with your system. See Section 10 System Care for replacement instructions.
- Incoming water supply line pressure to the system must be 20 psi (138 kPa) to 40 psi (275 kPa). If your dental water system’s supply line pressure is above 40 psi, install a water pressure regulator on the water supply line to your Cavitron Plus Ultrasonic Scaler.
- A manual shut-off valve on the dental water system supply line should be used so that the water can be completely shut-off when the office is unoccupied.
- In addition to the water filter supplied, it is recommended that a filter in the dental water system supply line be installed so that any particulates in the water supply will be trapped before reaching the Cavitron system.
- After the above installations are completed on the dental water supply system, the dental office water line should be thoroughly flushed prior to connection to the Cavitron system.

- Incoming water temperature to the Cavitron System should not exceed 25°C (77°F). If needed a device should be installed to maintain a temperature within this specification, or a Cavitron DualSelect Dispensing System attached to allow this system to be operated as a closed water system.

7.2 Electrical Requirements

- Incoming power to the system must be 100 volts AC to 240 volts AC, single phase 50/60 Hz capable of supplying 1.0 amps.
- The system power should be supplied through the AC power cord provided with your system.
- **WARNING:** To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

7.3 Unpacking the System

Carefully unpack your Cavitron Plus Ultrasonic Scaler and verify that all components and accessories are included:

1. Cavitron® Plus™ System with Handpiece Cable Assembly with swivel
2. Water Line Assembly (Blue) with Filter and Quick Disconnect
3. Additional Water Line Filter
4. Detachable AC Power Cord
5. Cavitron® Tap-On™ Wireless Foot Pedal
6. “AA” Batteries (4-Pack)
7. Auxiliary Cable for Tap-On™ Foot Pedal
8. Steri-Mate® Grip (not shown)
9. Steri-Mate® Detachable Sterilizable Handpiece
10. Cavitron® Ultrasonic Inserts (quantity optional)
11. Efficiency Indicator for Cavitron Inserts
12. Literature Packet

7.4 System Installation

- The Cavitron Plus System is designed to rest on a level surface. Be sure unit is stable and resting on four feet.
- The Cavitron Plus System should not be positioned such that access to the power input and AC power cord are limited.


- Placing unit in direct sunlight may discolor plastic housing.
- The system has been equipped with a Tap-On™ Wireless Foot Pedal which was factory synchronized to operate with the system's base unit. If your office has more than one Cavitron system with Tap-On™ Technology, it is recommended that you mark the Tap-On™ Foot Pedal and base unit for easy reference as to which Tap-On™ Foot Pedal operates with which base unit. Should resynchronization be necessary, follow the instructions in section 7.10.

7.5 Power Cord Connection


- Verify the Main Power ON/OFF switch, located at the center front underside of the System, is set to the OFF (O) position before proceeding.

- Insert the AC power cord into the power input on the back of the System.
- Insert the pronged plug into an AC wall outlet.

7.6 Water Supply Line Connection

- Grasp the Water Supply Line (blue hose) by the end opposite the quick-disconnect and insert it into the water inlet connector until fully seated.

Press ring to release water supply tube.

- Connect the quick-disconnect to the dental office water supply or a Cavitron DualSelect Dispensing System.
- Inspect all connections to make certain there are no leaks.
- To remove the water line from the Cavitron Plus System, turn off the dental office water supply. Disconnect the water supply line from the dental office water supply. If a quick-disconnect connector is attached to the end of the hose, relieve the water pressure by pressing the tip of the connector in an appropriate container and allow water to drain. To remove the hose from the system, push on the outer ring of the system's water inlet and gently pull out the water line.

7.7 Tap-On™ Foot Pedal Battery Installation/Replacement

Look for blinking communications light.

- Turn Tap-On™ Wireless Foot Pedal over and using a Philips screwdriver carefully remove battery cover screw and battery cover. If applicable, remove used batteries and install two new “AA” batteries as shown. Do not depress Tap-On™ Foot Pedal while installing batteries.
- The communication light will blink for approximately two seconds to indicate the Tap-On™ Foot Pedal's ability to communicate with the unit. If the light does not blink, check the batteries. If the batteries are good and the light doesn't blink, a communications error may exist. To re-establish communication with Tap-On™ Foot Pedal review Synchronization procedure, section 7.10.
- The remote frequency communication can be bypassed using the auxillary Tap-On™ Foot Pedal cable. Refer to Section 11.2 Technical Support and Repair for further action.
- Replace the battery cover and screw and hand tighten cover with Philips screwdriver.
- Remove batteries if Tap-On™ Foot Pedal is to be stored for an extended period of time.

7.8 Tap-On™ Foot Pedal Synchronization

The Tap-On™ Wireless Foot Pedal supplied with your system has been factory synchronized with the base unit. Should a replacement Tap-On™ Foot Pedal be necessary, synchronization will be required prior to system operation. Perform the following steps to synchronize the Tap-On™ Foot Pedal with the base unit.

1. Turn the MainPower switch located at the center

front underside of the system to the OFF (O) position.

2. Install a new set of “AA” batteries into the foot control (see section 7.9). Leave the battery cover of the Tap-On™ Foot Pedal open so the red push button is accessible.
3. Maintain a distance of no more than 10 feet (3 meters) between the base unit and Tap-On™ Foot Pedal during the synchronization process.
4. Remove any inserts from the handpiece and adjust the Power Level Control out of Rinse Mode. Turn the Main Power switch to the ON (I) position and wait for the Diagnostic Display graphics to light (refer to Section 8.2).
5. While all graphics are lit, press the Purge button, located on the Diagnostic Display.


The graphics will begin to blink in a sequential pattern, representing the synchronization mode. This mode will last 5 to 6 seconds.

6. During this mode, press the red button located in the battery compartment of the Tap-On™ Foot Pedal. This will complete the synchronization process.

7. Synchronization is successful when all graphic lights blink at the same time.
8. To verify proper communication, press the foot control to the Boost position (Tap-On™ Foot Pedal fully depressed – 2nd position) and ensure the Boost graphic on base unit illuminates.
9. Attach battery cover and tighten the screw.
10. In the event communication cannot be established, temporarily use the supplied Auxillary Tap-On™ Foot Pedal Cable to connect the Tap-On™ Foot Pedal directly to the unit.

SECTION 8: Cavitron Plus Ultrasonic Scaler Description

8.1 System Controls

Ultrasonic Power Level Control

Turn knob to select the ultrasonic power level for operation. Turning the knob clockwise increases the distance the insert tip moves (stroke) without changing the frequency; turning the knob counterclockwise decreases the distance the insert tip moves (stroke) without changing the frequency.

The Blue Zone is an extended low-power range for improved patient comfort when subgingivally scaling.

Rinse

Turn the ultrasonic power level control knob fully counterclockwise until a “click” is heard. Rinse mode is for use during an ultrasonic scaling procedure when lavage is desired with minimal cavitation.

Diagnostic Display


See section 8.2 for more details.

Handpiece Holder

Securely holds the system’s handpiece, or cable connector when handpiece is not installed.

Steri-Mate® Handpiece

Accepts all Cavitron® 30K™ Ultrasonic inserts
See section 8.3 for more details.

Blue Zone

Rinse

Main Power ON/OFF Switch

ON/OFF Switch located at the center front underside of the system.

Tap-On™ Wireless Foot Pedal

Eliminates the need to hold down or pump foot pedal.
See section 8.6 for more details.

8.2 Diagnostic Display Indicators and Control

Service Indicator

Illuminates when the system is not functioning properly. This display has three distinct modes.

- A fast blink (3 blinks per second) indicates an improper set-up.
- A slow blink (1 blink per second) means the system is operating out of factory specifications.
- A steady light indicates the system is overheating.

Refer to section 11.1 for Troubleshooting guidelines.

Power Indicator


Illuminates (3 sec. delay) when the Main Power ON/OFF Control Switch is ON ("I" position).

Low Battery Indicator

Illuminates when the Tap-On™ Foot Pedal battery power is approaching end of life. Replace batteries as instructed in section 7.9.

Boost Indicator

Illuminates when the Boost Mode has been activated with the Tap-On™ Foot Pedal.

Turbo Mode


When pressed, extra power is delivered to the system, up to 25%. This will increase the stroke of the tip. Turbo power remains on until the button is pressed again or unit is turned off. (Purple arrows illuminate when Turbo Mode is on).

Purge Control

Illuminates when the Purge function is activated. To activate Purge, remove insert from the handpiece, turn the handpiece lavage control to maximum water flow and press the Purge button. Water will purge through system for 2 minutes. To deactivate mode during the 2 minute cycle, press Purge button again or press Tap-On™ Foot Pedal.

The Purge Control is also used during the Tap-On™ Foot Pedal Synchronization process. See section 7.10 for more details.

8.3 Handpiece / Cable

Lavage Control

Turn the Lavage Control to select flow rate during system operation. Flow rate is based on a scale from 1 to 6. Turn clockwise toward 6 to increase flow at insert tip. Turn counter-clockwise toward 1 to decrease flow. The flow rate through the handpiece also determines the temperature of the lavage. Lower flow rates produce warmer lavage. Higher water flow rates produce cooler lavage.

If the handpiece becomes warm, increase the flow rate. With experience the Dental Healthcare Professional will be able to determine the best flow rate setting for optimum operating efficiency and patient comfort.

Swivel Feature


Reduces cable drag as handpiece rotates during procedures.

Steri-Mate® Grip accessory (not shown)

The Steri-Mate® Grip is designed to provide an ergonomic and comfortable grasp of the handpiece. The grip is sterilizable and is available in several different colors as an accessory for your Steri-Mate Handpiece. See installation instructions provided with the grip.

8.4 Caviron 30K Ultrasonic Inserts

The many styles of Caviron and Caviron Bellissima 30K Ultrasonic Inserts are easily interchangeable for various procedures and applications. See enclosed literature for specific information.

8.5 Tap-On™ Technology Wireless Foot Pedal Operation

Using Foot Pedal in Tap-On™ Mode

For scaling procedures, Tap-On™ Technology eliminates the need to hold the pedal down. Tapping the foot pedal once activates ultrasonic power or rinse mode for approximately 4 minutes. Tapping the foot pedal while in Tap-On™ mode disables the ultrasonic power and water flow. Boost is still available while scaling in Tap-On™ mode. To use boost, simply depress the foot pedal to the second position (all the way to the floor) to activate and hold as long as boost is desired. Release foot pedal to return to Tap-On™ mode.

TIPS:

- Tap-On™ Technology will not run water unless an insert is in the handpiece
- A sensor in the handpiece holder will prevent Tap-On™ Technology from operating when the handpiece is in the holder.
- If the foot pedal is not tapped quickly, it will function in a conventional manner.

How to Disable and Enable Tap-On™ Technology

Tap-On™ Technology feature can be disabled by simultaneously holding the Purge and Turbo buttons for a period of approximately 5 seconds. The two buttons will blink approximately 6 times. When the buttons are released, they will blink an additional 6 times to confirm Tap-On™ has been disabled.

Tap-On™ Technology feature can be enabled by simultaneously holding the Purge and Turbo buttons for a period of approximately 5 seconds. The two buttons will blink approximately 6 times to confirm Tap-On™ has been enabled.

Using Foot Pedal without Tap-On™ Mode

For scaling operation, the first position activates both the ultrasonic energy and lavage at the insert tip. The second position activates the Boost Mode. The Boost Mode (fully depressed Tap-On™ Foot Pedal) increases the ultrasonic power level for quick removal of tenacious deposits without adjusting the power level knob. To deactivate Boost Mode, release Tap-On™ Foot Pedal to first position.

8.6 Accessories and User Replaceable Parts

8.6.1 Accessories

1. AC Power Cord
2. Tap-On™ Technology Wireless Foot Pedal
3. Auxiliary Tap-On™ Foot Pedal Power Cable
4. Cavitron Steri- MATE® Sterilizable Handpiece
5. Cavitron 30K Ultrasonic Inserts

6. Cavitron DualSelect Dispensing system
7. Steri- MATE® Grip
8. Steri- MATE® Light Lighted Sleeve

8.6.2 User Replaceable Part Kits

1. Cavitron Insert Replacement O-ring Kits, 12/packs
Part Number 62351 (black) for plastic and soft grips
Part Number 62605 (green) for metal grips
2. Handpiece Cable O-Ring, Part Number 79357
3. Lavage (Water) Filter, 10/Pack, Part Number 90158

For detailed information, contact your local DENTSPLY Professional Representative or authorized DENTSPLY Professional Distributor.


SECTION 9: System Setup, Operation and Techniques for Use

9.1 Handpiece Setup

- This handpiece is sterilizable. Refer to Infection Control Information Booklet for sterilization instructions prior to using handpiece.
- Connect the Handpiece to the Cable Assembly by aligning the electrical connections. If Cable Assembly does not seat into the handpiece, gently rotate the handpiece until contacts align, then fully insert handpiece.
- Hold empty handpiece in a semi-upright position over a sink or drain. Activate the Tap-On™ Foot Pedal until water exits to release any air bubbles that might be trapped inside the handpiece. NOTE: Tap-On™ Technology only operates when an insert is in the handpiece.
- Lubricate the O-ring on the insert with water before placing it into the handpiece. Fully seat insert with a gentle push-twist motion. DO NOT FORCE.
- Turn the Lavage Control to select flow rate during system operation. Flow rate is placed on a scale from 1 to 6. Turn control clockwise toward 6 to increase flow at insert tip. Turn control counter-clockwise toward 1 to decrease flow. The flow rate through the handpiece also determines the temperature of the lavage. Lower water flow rates produce warmer lavage. Higher flow rates produce cooler lavage. If the handpiece becomes

warm, increase the flow rate. With experience the Dental Healthcare Professional will be able to determine the best flow rate setting for optimum operating efficiency and patient comfort.

9.2 Turbo Mode

Pressing the “Turbo” button on the scaling unit increases the unit’s ultrasonic power up to 25%. When you need more power for an extended period of time, simply press the Turbo button on the display panel (arrows will illuminate to show you are in “Turbo” mode).

DENTSPLY recommends that the clinician familiarize themselves with the available power levels throughout the power knob rotation in both normal and turbo modes. To do this, simply hold your favorite insert over the sink and adjust the power knob while observing the inserts spray pattern and toggling between both normal and turbo modes.

9.3 Boost Mode

Boost provides a temporary increase in ultrasonic scaling power for quick removal of tenacious calculus without touching the unit. Boost is activated by fully depressing the Tap-On™ Foot Pedal to second position (all the way to the floor). When boost is activated, the boost icon will illuminate on the display panel. Boost remains on as long as the clinician has the foot pedal pressed all the way down. In order to deactivate boost, release the Tap-On™ Foot Pedal to first position.

9.4 Patient Positioning

For optimal access to both the upper and lower arches, the backrest of the chair should be adjusted as for other dental procedures. This assures patient comfort and clinician visibility.

Have the patient turn his/her head to the right or left. Also position chin up or down depending upon the quadrant and surface being treated. Evacuate irrigant using either a saliva ejector or High Volume Evacuator (HVE).

9.5 Performing Ultrasonic Scaling Procedures

Note: Refer to the *Infection Control Information booklet supplied with your system and Section 10 of this manual for general procedures to be followed at the beginning of each day and between patients.*

- Follow precautions listed in the General and Ultrasonic sections of 4.2 Procedural Precautions
- The edges of Cavitron Ultrasonic Inserts are intentionally rounded so there is minimal danger of tissue laceration with proper ultrasonic scaling technique. Whenever the insert tip is placed in the patient’s mouth, the lips, cheek and tongue should be retracted to prevent accidental (prolonged) contact with the activated tip.
- Turn Power Level Control to select ultrasonic power level for operation. Clockwise increases system power. Power level will increase throughout the full range of the control. Hold the handpiece over a sink or drain. While in Tap-On™ mode, simply tap the Tap-On™ Foot Pedal to activate the system. (If Tap-On™ mode is turned off, press and hold the Tap-On™ Foot Pedal down to activate the system.) Check water spray to verify fluid is reaching the working end of the insert tip. Adjust the water lavage control until the water (lavage) flows with a rapid drip or small spray. Higher water flow settings provide cooler irrigation.
- It may be necessary to adjust lavage with the system in “Boost” mode (Tap-On™ Foot Pedal fully depressed) so adequate fluid will be available to cool tip to tooth interface.
- In general, it is suggested that a “feather-light-touch” be used for ultrasonic scaling. The motion of the activated tip and acoustic effects of the irrigating fluid, in most cases, are adequate to remove even the most tenacious calculus.
- Periodically check the Cavitron Ultrasonic Insert for wear with the Cavitron Insert Efficiency Indicator.
- The use of a saliva ejector or High Volume Evacuator (HVE) is recommended during all procedures.
- Set the system’s Power Level Control to the lowest efficient power setting for the application and the selected insert
- Keep the foot pedal near your foot to make it convenient to access.

9.6 Patient Comfort Considerations

Reasons for sensitivity

- Incorrect tip placement. The point should never be directed toward tooth root surfaces.
- Not keeping tip in motion on tooth. Do not allow the insert to remain in a static position on any one area of the tooth. Change the insert’s path of motion.
- Applying excessive pressure. Use a very light grasp and pressure, with a soft tissue fulcrum whenever possible, especially on exposed cementum.
- If sensitivity persists, decrease power setting and/or move from the sensitive tooth to another and then return.

SECTION 10: System Care

It is recommended that you perform the following maintenance procedures.

10.1 Daily Maintenance

START-UP PROCEDURES AT THE BEGINNING OF THE DAY:

1. Open the manual shut-off valve on the dental office water supply system.
2. Install a sterilized Steri-Mate® Handpiece onto the handpiece cable.
3. Set the Power Level Control to minimum and the Lavage Control to maximum.
4. Turn the system ON.
5. Hold the sterilized handpiece (without insert installed) over a sink or drain. Activate the Purge Control button.
 - The Purge button will light for two minutes indicating proper activation of the purge function.
 - If the Purge button is activated with an insert present in the handpiece, the button will blink for 3 seconds and disable. Remove the insert from the handpiece and press the Purge button again.
 - The Purge function can be interrupted at any time during the two minute cycle by pressing the Purge button again or by pressing the Tap-On™ Foot Pedal.
6. After completing the purge cycle, place a sterilized 30kHz Cavitron® Ultrasonic Insert into the handpiece and set the Power Level Control and Lavage Control to your preferred operating position for ultrasonic scaling.

BETWEEN PATIENTS:

1. Remove the used Cavitron® Ultrasonic Insert. Clean and sterilize following the Infection Control Procedures that were enclosed with your insert.
2. Hold the handpiece over a sink or drain and activate Purge function as described in Step 10 of the Start-Up procedure.
3. After the purge cycle is complete, turn the System to the OFF (O) position.
4. Remove the Steri-Mate® handpiece, clean and sterilize following the procedures outlined in the Cavitron Systems Infection Control Procedures booklet that was enclosed with your system.
5. Disinfect the surfaces of the cabinet, Power Cord, Handpiece Cable, Tap-On™ Foot Pedal and Cable assembly (if applicable), Water Supply and Air Supply lines as instructed in the Cavitron Infection Control Procedures booklet. **DO NOT SPRAY DISINFECTANT SOLUTION DIRECTLY ON SYSTEM SURFACES.**

6. Inspect the handpiece cable for any breaks or tears.
7. If using a closed water supply or DualSelect Dispensing System, check for adequate fluid volumes for the next patient.
8. When ready to use, place a sterilized Steri-Mate® Handpiece onto the handpiece cable assembly and insert a sterilized ultrasonic insert into the handpiece and adjust system controls as preferred.

SHUT-DOWN PROCEDURES AT THE END OF THE DAY:

Follow the “Between Patients” maintenance procedures, Steps 1 through 6. In addition, it is recommended to close the manual shut-off valve on the dental water supply system.

***NOTE: Broad spectrum hospital grade water-based disinfection solutions are preferred, such as Lysol IC. Some alcohol-based disinfectant solutions may be harmful and may discolor plastic materials.**

10.2 Weekly Maintenance

It is strongly recommended that this system be disinfected by chemically flushing the waterlines with a 1:10 Sodium Hypochlorite solution (NaOCl) at the end of each week. This can be accomplished by connecting this device to the Cavitron DualSelect Dispensing System or a number of other devices available from your local distributors. Where this device is connected to the Cavitron DualSelect Dispensing System, please follow the DualSelect system’s Directions for Use manual. If connected to another device, please follow those directions for use, keeping in mind that a chemical flush should be performed at maximum water flow for at least 30 seconds. The system should be left undisturbed for 10 minutes but no more than 30 minutes to allow the sodium hypochlorite solution to soak in the lines. As a suggestion, it is recommended that a sign be placed on the system stating that the SYSTEM IS BEING DISINFECTED WITH A STRONG DISINFECTANT AND SHOULD NOT BE USED. When ready, flush system with clean water for at least 30 seconds or until sodium hypochlorite odor disappears. **ALL CHEMICALS MUST BE FLUSHED FROM THE SYSTEM BEFORE IT IS READY FOR PATIENT USE.**

10.3 Monthly Maintenance

WATER LINE FILTER MAINTENANCE: When the water line filter becomes discolored, the filter should be replaced to prevent reduced water flow to the Cavitron Plus system. A 10-pack of replacement filters is available by ordering Part Number 90158 from your local DENTSPLY Professional distributor.

1. Verify that the system is turned OFF.
2. Disconnect the water supply hose from the dental office water supply. If a quick disconnect connector is attached to the end of the hose, relieve the water pressure by pressing the tip of the connector in an appropriate container and drain the water.
3. Grasp the fittings on either side of the filter disk and twist counterclockwise. Remove the filter section from either side of the water hose.
4. Install the replacement filter onto the water hose fittings. The filter should be positioned to match up with the correct hose fitting.
5. Hand tighten one hose fitting onto filter in a clockwise direction. Tighten second hose onto filter in clockwise direction. Reconnect the water supply line, operate unit to bleed the air and test for leaks.

SECTION 11: Troubleshooting

Although service and repair of the Cavitron Plus Ultrasonic Scaler should be performed by DENTSPLY personnel, the following are some basic troubleshooting procedures that will help avoid unnecessary service calls. Generally, check all lines and connections to and from the System. A loose plug or connection will often create problems. Check the settings on the System's controls.

11.1 Troubleshooting Guide

Symptom:

System operates: Tap-On™ Technology is not working

1. Tap-On™ Technology might be disabled. Refer to section 8.5.
2. Check to see if handpiece is in holder. Tap-On™ Technology is disabled when handpiece is in holder.
3. Check to see if the insert is secured inside the handpiece. Tap-On™ Technology is disabled when there is no insert in the handpiece.

Symptom:

System will not operate: No Power ON indicator

1. Check that the Main Power Switch is in the ON (I) position, and that the detachable Power Cord is fully seated in the receptacle on back of System.
2. Check that the system's power cord plug is fully seated in an approved AC wall outlet.
3. Check that the wall outlet is functional.

Symptom:

System will not operate: Power ON Indicator is illuminated

1. If the office has more than one Tap-On™ Foot Pedal, test each to ensure that the proper Tap-On™ Foot Pedal is being used. With a handpiece and insert installed, depress the Tap-On™ Foot Pedal to the first position. The system should dispense water. If none of the Tap-On™ Foot Pedals operate the system,

continue to the next step.

2. Resynchronize one Tap-On™ Foot Pedal to the system (see Section 7.8 Tap-On™ Foot Pedal Synchronization).

Symptom:

System operates: No water flow to insert tip or handpiece overheats

1. Assure that handpiece lavage control is properly adjusted.
2. Check for clogged insert. Replace insert if necessary.
3. Check that dental office water supply valves are open.
4. If the system is connected to DualSelect Dispensing System, check that fluid level in the selected bottle is sufficient. Make sure valves are open when using external water source.
5. Check that the water line filter is clean. Replace filter if needed.

Symptom:

System operates: No insert cavitation

1. Check that the Power Level Control is not in Rinse Mode.
2. Check the insert for damage and that it is properly installed in the handpiece.
3. Check that the handpiece is properly installed to the cable assembly.
4. Verify that the soft nozzle grip is flush with the hard plastic of the insert port.
5. Turn the system's Main Power Switch to the OFF (O) position. Wait 5 seconds and turn the system back ON.
6. If problem still exists, replace both "AA" batteries in Tap-On™ Foot Pedal with new "AA" batteries (Refer to section 7.7) or connect the Auxiliary Tap-On™ Foot Pedal Cable.

Symptom:

System operates: Purge Mode will not function – icon flashing

1. Check that there is no insert in the handpiece.
2. Check that handpiece is properly installed to the cable assembly.

Symptom:

System operates: Service Indicator blinking

- Fast Blinking (3 blinks per second)
 - Indicates improper set-up
1. If insert is in the handpiece, remove. Verify the handpiece is properly seated and hold the foot control for 2 seconds. If blinking stops, the system is ready for use. If blinking remains, continue to the next step.
 2. Attach a NEW handpiece and hold Tap-On™ Foot Pedal for 2 seconds. If blinking stops, the system is ready for use. Discard the old handpiece or return if within warranty. If blinking remains, continue to the next step.
 3. Install and fully seat an insert into handpiece. Hold Tap-On™ Foot Pedal for 2 seconds. If blinking stops, the system is ready for use. If blinking remains, continue to the next step.
 4. Install and fully seat a NEW insert in handpiece and

hold Tap-On™ Foot Pedal for 2 seconds. If blinking stops, system is ready for use. Discard old insert or return if within warranty. If blinking remains, refer to Section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.

- Slow Blinking (1 blink per second)
 - The system is not operating within factory specifications.
 1. Remove insert.
 2. Turn Main Power Switch OFF, (O) position. Wait five seconds. Turn unit ON, (I) position.
 3. Operate Purge function.
 4. If service indicator still blinks, refer to section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.

Symptom:

System operates: Service Indicator illuminated

1. Ensure that the base unit has adequate ventilation and is not near a heat source (i.e. radiator, heat lamp, sunlight or other heat producing operatory equipment).
2. Turn Main Power Switch to the OFF (O) position. Allow system to cool for 10 minutes and turn system to the ON (I) position. Verify light is not illuminated.
3. If light is still illuminated, refer to section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.

11.2 Technical Support and Repairs

For technical support and repair assistance call DENTSPLY Professional Cavitron CareSM Factory Certified Service at 1-800-989-8826 Monday through Friday, 8:00 A.M. to 5:00 P.M. (Eastern Time). For areas outside the U.S., contact your local DENTSPLY Professional representative.

**SECTION 12:
Warranty Period**

The Cavitron Plus Ultrasonic Scaler is warranted for TWO YEARS from date of purchase. The Steri-Mate® Handpiece enclosed with your system is warranted for SIX MONTHS from date of purchase. Refer to the Warranty Statement Sheet furnished with your system for full Warranty Statement and Terms.

**SECTION 13:
Specifications**

Electrical Voltage	Continuous (100-240 VAC)
Current	1.0 Amperes, Maximum
Phase	Single
Frequency	50/60 Hertz
Water Pressure	20 to 40 psig (138 to 275 kPa)
Water Temperature	< 25°C (77°F)
Water Flow Rate	Minimum Setting (CCW) < 15 ml/min Maximum Setting (CW) > 55 ml/min
Weight	4.4 lbs (2 Kg)
Dimensions	Height: 6 in. (15,24 cm) Width: 9.5 in. (24,13 cm) Depth: 8 in. (20,32 cm) Handpiece Cable length: 6.5 ft. (2.0 M) Auxillary Footswitch Cable length: 8 ft. (2.4 M) Water Supply Line length: 8 ft. (2.4 M)
Footswitch	Protection Class IPX1. Not for operating theatres.
Remote Communication	Frequency: 2405 to 2480 MHz Power: < 1mW Channels: 16
Operating Environment	Temperature: 15 to 40 Deg Celsius (59 to 104 Deg Fahrenheit) Relative Humidity: 30% to 75% (non-condensing)
Transport and Storage Conditions	Temperature: -40 to 70 Deg. Celsius (-40 to 158 Deg. Fahrenheit) Relative Humidity: 10% to 100% (non-condensing) Atmospheric Pressure: 500 to 1060 hPa

SECTION 14: Symbol Identification

AC POWER

TYPE B APPLIED PART EQUIPMENT

PROTECTIVE EARTH (GROUND)

IPX1

Footswitch not for operating theatres
Protection Class- IPX1
IPX1 Classification of ingress of water

Consult Instructions for Use

O/I

AC Power Switch (O = Off, I = On)

Footswitch

This symbol is a mandatory marking for devices entering the European market to indicate conformity with the essential health and safety requirements set out in European Directives. The symbol may be accompanied by a four-digit identification number of the notified body.

MEDICAL EQUIPMENT

WITH RESPECT TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL 60601-1/CAN/CSA-C22.2 No. 601.1, ANSI/AAMI ES60601-1 (2005, 3rd ed.) CAN/CSA-C22.2 No. 60601-1 (2008), 13VA

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

1) this device may not cause harmful interference, and
2) this device must accept any interference received, including interference that may cause undesired operation.

Foot pedal: FCC ID:TF3-DPD81861
IC: 4681B-DPD81861

Base: FCC ID:TF3-DPD81842
IC: 4681B-DPD81842

WEEE

Dispose of in accordance with the Waste Electrical and Electronic Equipment Directive 2002/96/EC of the European Parliament and the Council of the European Union

This is a wireless device.

SECTION 15: Classifications

- Type of protection against electric shock:
- Degree of protection against electric shock:
- Degree of protection against the harmful ingress of water:
- Mode of operation:
- Degree of safety of application in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide:
- According to medical device directive:
IIA (rule 9)

Class 1
Type B
Ordinary
Continuous

Equipment not suitable for use in the presence of flammable anaesthetic or oxygen.

SECTION 16: Disposal of Unit

U.S. - Dispose of the system components in accordance with state and local laws.

EU - Dispose of in accordance with the Waste Electrical and Electronic Equipment Directive 2002/96/EC of the European Parliament and the Council of the European Union.

SECTION 17: Electro-Magnetic Compatibility Precautions

Guidance And Manufacturer's Declaration - Electromagnetic Emissions

The Ultrasonic Scaler model G136 is intended for use in the electromagnetic environment specified below. The customer or the user of the Ultrasonic Scaler should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Ultrasonic Scaler uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The Ultrasonic Scaler is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	Not Applicable
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	The Ultrasonic Scaler is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance And Manufacturer's Declaration – Electromagnetic Immunity

The Model G136 is intended for use in the electromagnet environment specified below. The customer or the end user of the Model should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance And Manufacturer's Declaration - Electromagnetic Emissions

The Model G136 is intended for use in the electromagnetic environment specified below. The customer or the user of the Model G136 should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance Level	Electromagnetic environment - guidance
Radiated RF IEC 61000-4-3	3 V/m 26 MHz to 2.5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Model G136, including cables, than the Recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>$d = 1.7 \sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model G136 is used exceeds the applicable RF compliance level above, the Model G136 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Model G136.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between Portable and mobile RF communications equipment and the model @ 3Vrms

The model G136 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model G136 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model G136 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{v_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	-	0.12	0.23
0.1	-	0.34	0.74
1	-	1.7	2.3
10	-	3.7	7.4
100	-	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Recommended separation distance between Portable and mobile RF communications equipment and the model @ 10Vrms

The model G136 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model G136 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model G136 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{v_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	-	0.035	0.07
0.1	-	0.11	0.22
1	-	0.35	0.7
10	-	1.1	2.2
100	-	3.5	7

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

QUICK REFERENCE GUIDE: TROUBLESHOOTING

SYMPTOM	ACTION TAKEN
System will not operate: No Power ON Indicator	<ol style="list-style-type: none"> 1. Check that the Main Power Switch is in the ON (I) position, and that the detachable power cord is fully seated in the receptacle on back of system. 2. Check that the system's power cord plug is fully seated in an appropriate AC wall outlet. 3. Check that the wall outlet is functional.
System will not operate: Power ON Indicator is illuminated	<ol style="list-style-type: none"> 1. If the office has more than one Tap-On™ Foot Pedal, test each to ensure that the proper Tap-On™ Foot Pedal is being used. With a handpiece and insert installed, depress the Tap-On™ Foot Pedal to the first position. The system should dispense water. If none of the Tap-On™ Foot Pedals operate the system, continue to the next step. 2. Resynchronize one Tap-On™ Foot Pedal to the system (see Section 7.10 Tap-On™ Foot Pedal Synchronization).
System operates: No water flow to insert tip or handpiece overheats	<ol style="list-style-type: none"> 1. Assure that handpiece lavage control is properly adjusted. 2. Check for clogged insert. Replace insert if necessary. 3. Check that dental office water supply valves are open. 4. If the system is connected to DualSelect Dispensing System, check that fluid level in the selected bottle is sufficient. Make sure valves are open when using external water source. 5. Check that the water line filter is clean. Replace filter, if needed.
System operates: No insert cavitation	<ol style="list-style-type: none"> 1. Check that the Power Control is not in Rinse Mode. 2. Check the insert for damage and that it is properly installed in the Handpiece. 3. Check that the handpiece is properly installed to the cable assembly. 4. Verify that the soft nozzle grip is flush with the hard plastic of the insert port. 5. Turn the system's Main Power Switch OFF, (O) position. Wait 5 seconds and turn the system back ON. 6. If problem still exists, replace both "AA" batteries in Tap-On™ Foot Pedal with new "AA" batteries (Refer to Section 7.9) or connect Auxiliary Tap-On™ Foot Pedal Cable.
System operates: Service indicator blinking	<ol style="list-style-type: none"> 1. Fast blinking (3 blinks per second) – Indicates improper set-up <ol style="list-style-type: none"> A. If insert is in the handpiece, remove. Verify the handpiece is properly seated and hold the Tap-On™ Foot Pedal for 2 seconds. If blinking stops, the system is ready for use. If blinking remains, continue to the next step. B. Attach a NEW handpiece and depress Tap-On™ Foot Pedal for 2 seconds. If blinking stops, the system is ready for use. Discard the old handpiece or return if within warranty. If blinking remains, continue to the next step. C. Install and fully seat an insert into handpiece. Depress Tap-On™ Foot Pedal for 2 seconds. If blinking stops, the system is ready for operation. If blinking remains, continue to the next step. D. Install and fully seat a NEW insert in handpiece and depress Tap-On™ Foot Pedal for 2 seconds. If blinking stops, system is ready for use. Discard old insert or return if within warranty. If blinking remains, refer to Section 11.2 Technical Support and Repairs to have unit serviced as soon as possible. 2. Slow blinking (1 blink per second) - The system is not operating within factory specifications. <ol style="list-style-type: none"> A. Remove insert. B. Turn Main Power Switch OFF, (O) position. Wait five seconds. Turn unit ON, (I) position. C. Operate Purge function. D. If service indicator still blinks, refer to Section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.
System operates: Service indicator illuminated	<ol style="list-style-type: none"> 1. Ensure that the base unit has adequate ventilation and is not near a heat source (i.e. radiator, heat lamp, sunlight or other heat producing operatory equipment). 2. Turn Main Power Switch OFF (O) position. Allow system to cool for 10 minutes and turn system to the ON (I) position. Verify light is not illuminated. 3. If light is still illuminated, refer to Section 11.2 Technical Support and Repairs to have unit serviced as soon as possible.
System operates: Purge mode will not function - icon flashing	<ol style="list-style-type: none"> 1. Check that there is no insert in the handpiece. 2. Check that Handpiece is properly installed to the cable assembly.

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United States of America

DENTSPLY Professional
Technical Service and Repair Department
1301 Smile Way
York, PA 17404-1785
Phone: (800) 989-8826 or (717) 767-8502

Deutschland

DENTSPLY DeTrey GmbH
De-Trey-Strasse 1
78467 Konstanz
Germany
Phone: 7531 583 0

France

DENTSPLY DeTrey
17 Michael FARADAY
78380 Montigny Le Bretonneux
France
Phone: (1) 30 14 77 77

Australia

DENTSPLY (Australia) Pty. Ltd
11-21 Gilby Road
Mount Waverley, Victoria 3149
Australia
Phone: (61) 3-9538-8280

United Kingdom

DENTSPLY Ltd.
Building 3
The Heights
Weybridge
Surrey
KT13 0NY
United Kingdom
Phone: (0) 800 072 3313

Italia

DENTSPLY DeTrey Italia s.r.l.
Via A. Cavaglieri, 26
I-00173 Roma
Italia
Phone: (06) 723 3626

Canada

DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3 Canada
Phone: (905) 851-6060

DENTSPLY
PROFESSIONAL

Manufactured by:
DENTSPLY Professional
1301 Smile Way
York, PA 17404-1785 USA

DENTSPLY DeTrey GmbH
De-Trey-Str. 1
78467 Konstanz
Germany

Imported and
Distributed by:
DENTSPLY Canada
Woodbridge, Ontario
L4L 4A3