
---

# PINDEX<sup>®</sup> SYSTEM

---

<b>EN</b>	Owner's Guide .....	3
<b>DE</b>	Benutzerhandbuch .....	8
<b>FR</b>	Guide d'utilisation.....	13
<b>ES</b>	Guia de uso .....	18
<b>IT</b>	Guida per l'utente .....	23
<b>NL</b>	Gebruiksaanwijzing .....	28
<b>SV</b>	Bruksanvisning .....	33
<b>DA</b>	Køber vejledning .....	37
<b>FI</b>	Käyttöopas.....	42
<b>NO</b>	Bruksanvisninger .....	46
<b>PT</b>	Guia do Proprietário .....	51
<b>EL</b>	Οδηγός Ιδιοκτήτη .....	56

## Pindex® System

Packaging includes:

Pindex Laser: Owner's Guide

Contents:

Shaft Lock Tool

Chuck Tool

Allen Wrench

Reamer

Flexible Rubber Mold Set:

Full Arch

Left Quadrant

Right Quadrant

Mounting Template

50 Dual Pins and Sleeves

100 Self-Articulating Pins

100 Self-Articulating Sleeves

50 Articulator Pins

Practice Kit:

50 Pin Sleeves, White

50 Pin Sleeves, Grey

50 Pins, Long

50 Pins, Short

Unit Components:

1. Drill motor/housing
2. Work table
3. Lens/laser holder
4. Work table lock
5. Shaft lock hole
6. Fuses
7. On/off switch
8. Power cord connector
9. Lamp access cover slot
10. Red lens
11. Hand grip ring
12. Light beam set screw
13. Light beam switch
14. Knurled drill chuck
15. Carbide drill
16. Dust ring
17. Drill depth adjustment ring
18. Lamp Connector

Thank you for choosing Coltène/Whaledent's Pindex® pin-setting system. It has been designed and manufactured to provide you with many years of reliable service.

The Pindex System is designed to assist in fabricating accurate models and dies for use in the fabrication of dental restorations.

The Pindex System is used to accurately drill parallel holes into the underside of the model. Precision Pindex Pins may then be cemented into the parallel holes. Pins are selected based on space available.

Precision pinning of dies and sections of the master model allows the technician to disassemble and reassemble the model repeatedly and accurately to facilitate the fabrication of restorations.

This owner's guide will provide step-by-step instructions for the production of accurate models and dies.

### I. Pre-drilling instructions

1. Master model is poured in stone or epoxy (Fig. 1).
2. Model base is trimmed perfectly flat on model trimmer (Fig. 2), so that the model thickness from base to preparation margins is approximately 15 mm (Fig. 3).
3. Sections and dies are pre-planned and pin locations are marked with a pencil on tooth side of model (Fig. 4) prior to drilling holes (Fig. 5).
4. When two pins are in close proximity, it is important to allow sufficient space for the placement of the sleeves (Fig. 6). If this is not possible, a dual pin or Tri-Plus pin should be used (Fig. 7).
5. Each die or section must have either two pins or a dual pin (Fig. 13).

This product has been manufactured with an IEC 320 power cord inlet connector using a detachable cordset (line cord) to plug in at the rear of the unit. Please be sure to plug in the cordset into the inlet connector before the next step (Fig. 45).

### II. Drilling pin holes

1. For best results use a slightly moist model.
2. The power cord must be plugged into a properly grounded outlet. The power switch is located in the rear alongside the IEC 320 power cord inlet connector. The "on" position is designated with an "I". The "off" position is designated with an "O" (Fig. 46). Turn switch on; red pilot light will glow brightly.
3. Place model on work surface and align pencil marks with the illuminated dot (Fig. 5).

Firmly hold the model with both thumbs and grasp the large black ring with the remaining fingers. Slowly squeeze with even pressure, allowing entire work table housing to slide downward until the end stop is reached (the light will automatically switch off). Sequentially drill all holes to full depth. Do not force; advance slowly and let the drill do the cutting.

4. After all holes have been drilled, each hole which is to receive a dual pin or Tri-Plus is notched lingually with a thin separating disc (Fig. 8, 9 and 10).
5. The model and holes are thoroughly cleaned of debris with compressed air.

## II. Placing pins and sleeves

1. Place a small amount of cyanoacrylate cement into a depression in a small piece of soft wax (Fig. 11).
2. Each Pindex Pin is lightly dipped into the cement and seated into the appropriate pin hole. For convenience and accessibility, the short Pindex Pins and dual pins are cemented first (Fig. 12). The dual pin must be aligned so that the anti-rotational rib fits into the groove in the model (Fig. 12). The long Pindex Pins are then cemented into place (Fig. 13).
3. The grey sleeves are placed on the short pins (Fig. 14).
4. The metal sleeve is placed on the dual pin (Fig. 15).
5. The white sleeves are placed on the long pins (Fig. 16). When used, the brown or yellow sleeves are placed over the medium pins.
6. The flat part of the sleeve should face the flat part of the other sleeve in each die.

## IV. Pindex die model base fabrication

1. Place carding wax over the extension of the long pins. It should be wide enough so that it can be positioned to close over the sleeve openings of the short pins and the dual pins (Fig. 17 and 18).
2. Insert mounting template pins through holes in rubber mold. Place an articulating sleeve onto each pin.
3. a. For patient and/or doctor identification, insert aluminum I.D. strip provided (Fig. 19) into the grooves with typed side facing the mold wall and the extensions bent inward (Fig. 20). The mold is now ready for base fabrication (Fig. 21).
  - b. Laboratory name and logo can be displayed by permanently cementing a rubber stamp into the recess provided in the flexible mold.
4. Lubricate model.
5. Carefully vibrate stone under carding wax and around pin sleeves (Fig. 22).
6. Fill rubber mold with stone (Fig. 23).
7. Carefully seat die model to place (Fig. 24) and wipe away excess stone. Allow to set.
8. Separate articulating base plate from rubber mold (Fig. 25) and separate model from rubber mold.
9. Trim excess stone on model trimmer (Fig. 26).
10. Snap articulating pins into sleeves (Fig. 27).
11. Lubricate and mount on articulator (Fig. 28).
12. Separate model from articulator by pulling up firmly (Fig. 29).
13. Strip away carding wax (Fig. 30) so that all pins and/or sleeves are exposed (Fig. 31).
14. This creates a window which provides access to the Pindex Pins, so that the dies can be pushed up for easy removal from model (Fig. 32).
15. The finished master model with the pins can be removed in one piece by pushing up on the pins using the corner of the lab bench (Fig. 33).

16. Dies can be sectioned from the underside of the model using the Pindex Handsaw (Fig. 34) and then snapping the final 2 mm of stone, preserving the interproximal margins (Fig. 35).
17. Dies can be sectioned in the conventional manner from above using a fine blade in the Pindex Handsaw.
18. Dies are then individually trimmed and reassembled, ready for wax pattern fabrication (Fig. 36).

## V. Troubleshooting, maintenance and service

Proper maintenance and care of your Pindex System drilling unit will eliminate most of the cause for wear, replacement of parts or the need for service calls.

### Oiling and cleaning

All bearings throughout the machine are permanently lubricated and require no additional servicing.

Frequently vacuum or brush away all accumulated debris and dust. **DO NOT USE AIR PRESSURE TO CLEAN THE MACHINE.** The pressure may force debris into moving parts.

The accumulation of grit and debris and infrequent cleaning are the major causes for downtime, parts wear and frequent service calls.

Debris created during drilling operations must be frequently removed. Brushing or vacuuming will adequately clear away all debris, especially from the work table and drill chuck assembly area.

Surrounding the drill chuck assembly is a dust ring collector (Fig. 37). This ring collects a majority of the drilling debris and prevents its migration to other parts of the apparatus. Carefully vacuum or brush away the debris accumulation.

A clean instrument is the single most important assurance of continued accuracy.

### Excessive pressure required

If a great deal of pressure is required to hold down the model during drilling, inspect the drill carefully. It may be worn and require replacement.

Do not lower the work table suddenly or too quickly. Use a slow, gradual motion and always let the drill do the cutting.

### Drill replacement

1. Disconnect plug from wall power outlet.
2. Depress work table lock and gently elevate until free from drill/motor housing (Fig. 38, arrow 2).
3. Disconnect lamp connector leads.
4. Place shaft lock rod (Fig. 39) into shaft lock hole (Fig. 40, arrow 4) as far as it will easily go. Gently rotate until the shaft lock rod slides into the hole provided. This operation will lock the drill chuck in position and prevent rotation.
5. Place chuck rod into the chuck lock hole (Fig. 37, arrow 5).
6. To open chuck, rotate chuck rod and knurled drill chuck counter-clockwise until resistance is felt. Maintain pressure while removing drill (Fig. 42).
7. Remove drill by hand. Replace with new drill, pressing down firmly with a piece of soft wood until fully seated.
8. Rotate chuck rod and knurled drill chuck clockwise to lock drill firmly in position.

- Remove both locking rods from their respective holes, connect lamp leads and replace work table housing onto motor housing. Plug instrument into wall receptacle and proceed with work. Carefully store both locking rods in a safe location for future use.

#### Drilling depth adjustment

Drill depth has been adjusted at the factory. On rare occasions, after changing drills, it may be necessary to make a drill depth adjustment. This knurled, internally threaded hand nut may be raised or lowered until the proper drilling depth has been achieved (Fig. 43).

#### Light beam spotter adjustment

- The light beam spotter switch is activated by the hex socket set screw mounted in the work table surface (Fig. 44, arrow 12).
- With the hex key provided, raise or lower the set screw until the light beam spotter goes off at the moment that the full drilled depth is achieved.

#### Unit does not drill parallel holes

This is usually caused by:

- Debris accumulation on the work table.
- A model that is not flat along its entire bottom work surface.
- Allowing the model to rise up off the work surface against the drilling pressure. Hold the model down against the work table more firmly and raise drill more slowly.
- The model is rocked during drilling. This will cause drill breakage.
- The model is permitted to slide during drilling. This will also cause drill breakage.

#### Does not start

- Check to determine that plug has not slipped out of wall receptacle and that the power switch breaker is ON.
- The motor thermal overload protector may be activated. Allow 1/2 hour to cool.
- The fuse may have blown.

#### Fuse blows

- Turn the main power switch to OFF position.
- Remove the power cord from the wall outlet and the power cord receptacle.
- The fuse compartment is located adjacent to the power cord receptacle. To access, lightly insert a small, flat-headed screwdriver in the slot on the upper section of the fuse compartment. Apply slight outward pressure to slide out the fuse compartment.
- Gently lift the center tab of the fuse compartment to release the clear plastic fuse holder. Slide the fuse holder completely out of the fuse compartment (Fig. 47).
- Examine the fuses. If either or both are blown, replace with fuses that match the unit's voltage. If both fuses appear to be in good condition and power is not activated when the main power on/off switch is turned ON, return the unit to your dealer or Coltène/Whaledent. See section IX Unit Specifications.
- To re-assemble the fuse compartment, slide the fuse holder back into the fuse compartment until it snaps into place. Replace the fuse compartment into the lower portion of the

power receptacle; it will snap into place.

- Make sure the main power switch is OFF before you plug the power cord into the power cord receptacle and the wall outlet.

#### Electric motor stops

The motor is provided with an automatic resetting inaccessible thermal overload protector. Always turn the unit OFF when not in actual use.

In the event the unit is accidentally left ON for a long period of time, the drill motor may automatically cut OFF.

- Switch the unit OFF.
- Allow to cool for approximately 1/2 hour.
- Turn ON and continue operation.

#### Squeal

The usual cause for squeal is grit and debris embedded in rotating parts. Turn OFF the machine, brush and vacuum away all debris accumulation in the rotating drill chuck area.

## VI. Accessory use directions

### A. Long, Medium, Short Pins and Sleeves

See Section III, Step 3 – Placing Pins and Sleeves

### B. Dual Pins and Sleeves

See Section III, Step 2 – Placing Pins and Sleeves

### C. Self-Articulating Pins and Sleeves

See Section IV, Step 2 – Pindex Die Model Base Fabrication

### D. Flexible Molds

See Section IV, Step 2 – Pindex Die Model Base Fabrication

### E. Pindex System Handsaw

See Section IV, Step 16 – Pindex Die Model Base Fabrication

### F. Pindex Laser Retro-fit Kit

#### Installation Instructions

- Tape a small piece of white paper to the work table to cover the drill hole. Do not remove paper until installation is complete.
- Turn the Pindex unit ON.
- Drill a hole through the paper.
- Turn the unit OFF after drilling the hole.
- Remove the work table from the motor housing by pressing in the work table lock-pin and lifting the work table away from the motor housing.
- Disconnect the black and white wires from the light arm to the motor housing by gently pulling them out of the connectors.
- Remove the light arm from the work table by removing the hex cap screws on the underside of the work table. Use the hex key provided.
- Place the new laser-arm into approximately the same position as the old light arm. You can use the same screws and washers or the screws and washers provided. Fasten the screws snug, not tight. The laser arm must be adjusted into a final position by hand.

9. Push the lamp wires gently into the connectors and turn the unit ON. It does not matter which wire goes into which socket.

**CAUTION:** Keep your hands away from the drill while the unit is on. Do not stare into the laser beam.

10. Move the laser-arm so that laser beam is pointing into the hole in the white paper on the work table.
11. Turn the unit OFF and tighten the screws underneath the platform to secure the arm into position.
12. Re-align the work table housing and the motor housing and push in the work table lock-pin. Slowly push the platform down until the work table lock-pin clicks into place. Turn the unit ON and check if the laser beam aligns with the hole in the paper. If the laser beam does not align with the hole, repeat steps 8 through 11.
13. Remove the paper from the work table. Adhere the two labels provided to the back of the unit, directly above the serial number label. The laser light should illuminate the tip of the drill when turned on.

## VII. Accessories

PX102	100	Medium Pins and Sleeves
PX116	100	Long Pins and Sleeves
PX117	100	Short Pins and Sleeves
PX114S	1000	Long Pins and Sleeves
PX112S	1000	Medium Pins and Sleeves
PX115	1000	Short Pins and Sleeves
PX145	100	Self-Articulating Pins and Sleeves
PX144	1000	Self-Articulating Pins and Sleeves
PX145S	1000	Self-Articulating Sleeves
PX146	3	Flexible Molds
PX147	1	Flexible Mold – Full Arch
PX148	2	Flexible Molds – Left and Right Quadrant
PX200		Pindex System Handsaw and Blades
PX125		Pindex Laser Retro-fit Kit

## VIII. Replacement parts

PX119	3	Carbide Drills
PX120	1	Reamer
PX124	10	Bulbs
PX201	12	Replacement Blades for Handsaw, 0.007" / 0.18mm
PX202	12	Replacement Blades for Handsaw, 0.010" / 0.25mm
PX301	100	Replacement Blades for Handsaw, 0.007" / 0.18mm
PX302	100	Replacement Blades for Handsaw, 0.010" / 0.25mm
P38	12	Proxima® Saw Blades

## IX. Unit specifications

Environmental and electrical conditions

- Indoor use only
- Transient over voltages category II.
- Pollution degree 2.
- Class 1 equipment
- Supply voltage fluctuations not to exceed  $\pm 10\%$  of nominal

Transportation and storage conditions

- Temperature  $-40^{\circ}\text{C}$  to  $+70^{\circ}\text{C}$  ( $-40^{\circ}\text{F}$  to  $158^{\circ}\text{F}$ )

## Pindex Laser

Cat. No.	Net Weight	Voltage $\pm 10\%$	Current	Watts	Fuse
PX113L	8.25 lbs./ 3.7 kg.	115 60 HZ	0.7 amp nominal	100 nominal	2 amp slow blow
PXM113L	8.25 lbs./ 3.7 kg.	100 50/60 HZ	0.7 amp nominal	100 nominal	2 amp slow blow
P113LCE	8.25 lbs./ 3.7 kg.	220/230 50/60 HZ	0.5 amp nominal	100 nominal	3AG .75 amp
P113LUK	8.25 lbs./ 3.7 kg.	220/230 50/60 HZ	0.5 amp nominal	100 nominal	3AG .75 amp

Relative humidity 10% to 100% including condensation

- Atmospheric Pressure 500 hPa to 1060 hPa  
(7.25 psi to 15.4 psi)

Operational ambient conditions

- Temperature  $5^{\circ}\text{C}$  to  $+40^{\circ}\text{C}$  ( $41^{\circ}\text{F}$  to  $104^{\circ}\text{F}$ )
- Relative humidity 50% to 80%
- Atmospheric Pressure 700 hPa to 1060 hPa  
(10.2 psi to 15.4 psi)
- Altitude  $\leq 2000\text{m}$  (6561 ft)

## X. Product markings

Caution (refer to accompanying documents)

Caution, risk of electric shock


Protective Conductor Terminal

## XI. Cautions and Warnings

Note: Failure to comply with these Caution and Warning statements may result in injury to person or equipment damage.

 <b>WARNING!</b>	Use only supplied power cord or equivalent.
	Do not disassemble the unit beyond removal of the table. There are no user serviceable parts.
	Ensure fingers, loose clothing, and jewelry are clear of the drill before pressing the base down.
	Do not use air pressure to clean the machine.
	Unit should be turned off when not in use.
	Do not stare into laser beam.
	This product contains a laser system and is classified as a CLASS I LASER PRODUCT. This product complies with IEC 60825-1:2001, 21 CFR 1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007. To use this model properly, read the instruction manual carefully and keep this manual for your future reference. In case of any trouble with this model, please contact your nearest "AUTHORIZED service distributor." To prevent direct exposure to laser beam, do not try to open enclosure.

 <b>CAUTION!</b>	Do not subject unit to shock or impact.
	This product is intended for use by trained healthcare professionals.
	Use this product only in accordance with the instructions in the owner's guide.
	Unspecified or improper use of this product may impair safety protection.
	This equipment may cause radio interference or may disrupt nearby equipment. It may be necessary to take mitigation measures such as relocating the unit or shielding the location.

## XII. Warranty

Our products are carefully manufactured to meet stringent quality assurance requirements. Our products are manufactured from new parts or new and serviceable used parts. Regardless, our warranty terms apply. This product has been developed specifically for use in dentistry and is intended to be operated only by qualified dental professionals in accordance with the instructions contained in this guide. However, notwithstanding anything contained herein to the contrary, the user shall at all times be solely responsible for determining the suitability of the product for the intended purpose and the method of its use. Any guidance on application technology offered by or on behalf of the manufacturer, whether written, verbal or by demonstration, shall not relieve the dental professional of his/her obligation to control the product and to make all professional judgments regarding its use.

Our products are warranted in accordance with the terms of a written Certificate of Limited Warranty accompanying each product. Except for the warranties specifically set forth in the Certificate of

Limited Warranty, Coltène/Whaledent Inc. provides no warranties or guarantees of any kind covering the product, expressed or implied, including, without limitation, any warranties as to merchantability or fitness for a particular purpose. The purchaser/user is referred to the Certificate of Limited Warranty for all of the terms, conditions and limitations of the warranty covering this product. This Section of the user manual is not intended to in any way modify or add to the warranty provided in the Certificate of Limited Warranty.

Any claim for damage or breakage to the product in transit should be made to the carrier promptly upon discovery. C/W does not warrant the product against shipping damage.

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**Coltène/Whaledent GmbH + Co. KG**
**Coltène/Whaledent AG**