Better outcomes start here

Palodent®V3 Sectional Matrix System



Palodent[®] V3

"70% of clinicians find contact creation to be the most challenging part of a Class II restoration."

¹ DentalTown (2012). Restorative Dentistry. Monthly Poll: What is the most challenging part of Class II Restoration?




Producing high quality Class II Restorations

The Palodent^{*}V3 system is an innovative sectional matrix system with Nickel-Titanium rings, anatomically shaped matrix bands and adaptive self-guiding wedges that provides predictable, tight contacts and restorations that accurately replicate the natural tooth anatomy.

"Love it, love it. I can't believe I practiced for the last 13 years without it. I am embarrassed that my Class II fillings would look like straight walls without the contours!" M. Rodrigue, DDS The Palodent^{*}V3 matrix system delivers improved outcomes, even when faced with wide cavities and missing cusps. The contoured design of the matrices and wedges ensure there is minimal finishing of the restoration. Used by tens of thousands of clinicians in dental offices and schools globally, better dentistry starts with Palodent^{*}V3



Better outcomes with a Sectional Matrix System



Circumferential Band (Tofflemire)


- × Fails to restore proximal anatomy
- Thin contact at the marginal ridge
- Increased likelihood of fracture, occlusal interference, recurrent caries and periodontal disease
- Large food trap below
- Increased chance of light or open contact
- Extensive finishing may be required

- Sectional Matrix System Restoration
- Easy to use system
- Natural contours
- Tight, anatomically correct contact point
- Minimal interproximal finishing

Sectional System vs Circumferential System

The Tofflemire circumferential matrix band was developed in the 1940s to aid placement of dental amalgam. It was not designed with today's modern restorative materials in mind, which can result in less than optimal outcomes when placing composites.

In contrast, sectional matrix systems were developed specifically for placing composite materials.



Take the stress out of Class II restorations


- 1 Strong Ni-Ti Rings for outstanding spring strength and memory
- 2 Glass-fiber reinforced plastic tines are V-shaped to accommodate the wedge
- **3** Pin tweezer hole for easy placement & removal of wedge
- 4 Top tab for easy insertion and side tab for removal of matrix
- 5 Gingival apron to prevent gaps in gingival-axial corner
- 6 Greater curvature: matrix wraps around tooth
- 7 Pronounced marginal ridge for natural-looking anatomy




4 steps to better Class II's*




Using the Pin-Tweezers, place a Palodent[®]V3 Wedge or WedgeGuard firmly from the lingual or buccal side and proceed with prep.



Place Palodent[®]V3 Matrix with Pin Tweezers.



After prep, if using Palodent[®]WedgeGuard, support Wedge with mirror or dental instrument while using Pin Tweezers to detach the shield.



Place the Palodent[®]V3 Ring on top of the wedge to provide separation.

Getting started

Every component in the Palodent[®]V3 Sectional Matrix System is optimized to make every step easier. Request a demo to see first hand how the Palodent[®]V3 components stack up against your current system.



659700v

Introductory Kit

100 Matrices – 25 each size: 3.5 mm, 4.5 mm, 5.5 mm, 6.5 mm 75 Wedges – 25 each size: Small, Medium, Large 30 WedgeGuards – 10 each size: Small, Medium, Large 1 Universal Ring; 1 Narrow Ring; 1 Forceps; 1 Pin Tweezers

Universal Ring (light blue)

- Nickel-titanium for strength and memory
- Stackable for multiple restorations
- V-shaped, glass-fiber reinforced plastic tines

Narrow Ring (dark blue)

- Tines set closer together for optimal separation on smaller teeth
- Up to 4kgs of separating force even on premolars

EZ-Coat Matrices (5 sizes)

- Ultra-thin, non-stick coating makes matrix removal easier and reduces light reflection, improving depth of field for better visibility
- Greater curvature for anatomical precision and less interference with ring placement
- Sizes color-coded for easy identification

Matrices (5 sizes)

- Pronounced marginal ridge for accurate anatomy
- Gingival apron is optimally shaped to prevent gaps in gingival-axial corner

Wedges (3 sizes)

- Compresses on entry and flares upon exit for easy placement and seal
- Hollow underside allows placement of second wedge from opposite side and protection of the papilla

Wedgeguards (3 sizes)

- Guard detaches after prep, leaving the wedge in place
- Facilitates faster and easier cutting of cavity and lowers risk of interproximal damage to adjacent tooth

Pin Tweezers

- Pin-in tip positively grips hole in matrix and wedge for control
- Naturally closed position minimizes risk of dropping matrix band or wedge
- Built-in ball burnisher tip

REORDER INFORMATION

Ring Refills 659900v Universal (1)

659760v	Universal (2)
659770v	Narrow (2)

Matrix Refills EZ Coat

659610v	Size 3.5 mm (50)
659620v	Size 4.5 mm (50)
659630v	Size 5.5 mm (50)
659640v	Size 5.5 mm (90)
659650v	Size 6.5 mm (50)
659660v	Size 7.5 mm (50)

Matrix Refills

659710vSize 3.5 mm (50)659720vSize 4.5 mm (50)659730vSize 5.5 mm (50)659740vSize 5.5 mm (100)659750vSize 6.5 mm (50)659860vSize 7.5 mm (50)

Wedge Refills

659780v	Small (100)
659790v	Medium (100)
659800v	Large (100)

 WedgeGuard Refills

 659830v
 Small (50)

 659840v
 Medium (50)

 659850v
 Large (50)

Instrument Refills		
659810v	Forceps Refill (1)	
659820v	Pin Tweezers Refill (1)	









THE DENTAL SOLUTIONS COMPANY™



K79200238-00 © Dentsply Sirona 12/2016