

Trade name: Steribim® plus (REF 54923)

Product no.: 54923

Current version: 1.0.0, issued: 17.01.2022 Replaced version: -, issued: - Region: GB

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Trade name

Steribim® plus (REF 54923)

1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses of the substance or mixture

Manufacturing of dental prosthesis in a dental laboratory

Uses advised against

No data available.

1.3 Details of the supplier of the safety data sheet

Address

BEGO Bremer Goldschlägerei Wilh. Herbst GmbH & Co. KG Wilhelm-Herbst-Str. 1

28359 Bremen

Telephone no. +49/ 421/ 2028 - 0 Fax no. +49/ 421/ 2028 - 115 e-mail msds@bego.com

Information provided by / telephone

Research & Development Department - Materials, alloys and ceramics; +49/ 421/ 2028 - 130 (Chief Development Officer alloys)

Advice on Safety Data Sheet

msds@bego.com

1.4 Emergency telephone number

For medical advice (in German and English): +49 (0)551 192 40 (Giftinformationszentrum Nord)

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification information

This product does not meet the classification and labelling criteria given in the Regulation (EC) No 1272/2008 (CLP).

2.2 Label elements

Not relevant

2.3 Other hazards

PBT assessment

The components of this product are not considered to be a PBT.

vPvB assessment

The components of this product are not considered to be a vPvB.

SECTION 3: Composition/information on ingredients

3.1 Substances

Not applicable. The product is not a substance.

3.2 Mixtures

Hazardous ingredients

This product does not contain substances to be mentioned according to EU Regulation No 1907/2006 (REACH), annex II.

SECTION 4: First aid measures



Trade name: Steribim® plus (REF 54923)

Product no.: 54923

Current version: 1.0.0, issued: 17.01.2022 Replaced version: -, issued: - Region: GB

4.1 Description of first aid measures

General information

If the patient is likely to become unconscious, place and transport in stable sideways position. Remove contaminated clothing and shoes immediately, and launder thoroughly before reusing.

After inhalation

Remove affected person from the immediate area. Ensure supply of fresh air. In case of persisting adverse effects consult a physician.

After skin contact

When in contact with the skin, clean with soap and water. Get medical attention if pain still persists.

After eve contact

Remove contact lenses, irrigate copiously with clean, fresh water for at least 15 minutes, holding the eyelids apart and seek medical advice. Get medical attention if pain still persists.

After indestion

Rinse the mouth thoroughly with water. Never give anything by mouth to an unconscious person. Let plenty of water be drunk in small gulps. Call a doctor.

4.2 Most important symptoms and effects, both acute and delayed

No data available.

4.3 Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media

Carbon dioxide; Extinguishing powder; Water spray jet; Foam

Unsuitable extinguishing media

High power water jet

5.2 Special hazards arising from the substance or mixture

The product is not flammable.

5.3 Advice for firefighters

Use self-contained breathing apparatus. Wear full protective suit. Adapt extinguisher and fire-fighting measures to fire in the environment. Cool endangered containers with water spray jet. Fire residues and contaminated firefighting water must be disposed of in accordance with the local regulations.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

For non-emergency personnel

Refer to protective measures listed in sections 7 and 8. Ensure adequate ventilation.

For emergency responders

No data available. Personal protective equipment (PPE) - see Section 8.

6.2 Environmental precautions

Do not discharge into the drains/surface waters/groundwater. Do not discharge into the subsoil/soil.

6.3 Methods and material for containment and cleaning up

Pick up with absorbent material (e.g., sand, kieselguhr, acid binder, universal binder, sawdust). When collected, handle material as described under the section heading "Disposal considerations".

6.4 Reference to other sections

No data available.

SECTION 7: Handling and storage

7.1 Precautions for safe handling



Trade name: Steribim® plus (REF 54923)

Product no.: 54923

Current version: 1.0.0, issued: 17.01.2022 Replaced version: -, issued: - Region: GB

Advice on safe handling

Ensure adequate ventilation. Avoid formation of aerosols.

General protective and hygiene measures

Wash hands before breaks and after work. Do not eat, drink or smoke during work time. Remove soiled or soaked clothing immediately. Keep away from foodstuffs and beverages. Provide eye wash fountain in work area. Have emergency shower available.

Advice on protection against fire and explosion

Pay attention to general rules of internal fire prevention.

7.2 Conditions for safe storage, including any incompatibilities

Technical measures and storage conditions

Keep container tightly closed and dry in a cool, well-ventilated place.

Requirements for storage rooms and vessels

Keep only in the original container. Containers which are opened must be carefully resealed and kept upright to prevent leakage.

Incompatible products

Substances to be avoided, see section 10.

7.3 Specific end use(s)

No data available.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

No parameters available for monitoring.

8.2 Exposure controls

Appropriate engineering controls

Provide adequate ventilation. This should be achieved by the use of local exhaust ventilation and good general extraction.

Personal protective equipment

Respiratory protection

Breathing apparatus in the event of aerosol or mist formation. Particle filter half mask, filter P1 - P3

Eye / face protection

Safety glasses with side protection shield (EN 166)

Hand protection

Sufficient protection is given wearing suitable protective gloves checked according to i.e. EN 374, in the event of risk of skin contact with the product. Before use, the protective gloves should be tested in any case for its specific workstation suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn. Design operations thus to avoid permanent use of protective gloves.

Appropriate Material nitrile rubber

Material thickness 0.4 mm

Appropriate Material natural latex

Material thickness 0.4 mm

Other

Normal chemical work clothing.

Environmental exposure controls

Avoid release into sewage and environment.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

State of aggregation
quid



Trade name: Steribim® plus (REF 54923)

Product no.: 54923

Current version: 1.0.0, issued: 17.01.2022 Replaced version: -, issued: - Region: GB

Form/Colour
paste
colourless; whitish

Odour
nearly odourless

pH valueValue< 8</td>Sourcesupplier

Boiling point / boiling range

No data available

Melting point/freezing point

No data available

Decomposition temperature

No data available

No data available

Flash point

Ignition temperature

No data available

Oxidising properties

Not applicable

Flammability
Not applicable

Source supplier

Lower explosion limit

No data available

Upper explosion limit

No data available

Vapour pressure

No data available

Relative vapour density
No data available

140 data available

Relative density
No data available

Density

No data available

Solubility

No data available

Partition coefficient n-octanol/water (log value)

No data available

Viscosity

No data available

Particle characteristics

No data available

9.2 Other information

Other information

No data available.

SECTION 10: Stability and reactivity



Trade name: Steribim® plus (REF 54923)

Product no.: 54923

Current version: 1.0.0, issued: 17.01.2022 Replaced version: -, issued: -Region: GB

10.1 Reactivity

Dangerous reactions are not expected if the product is handled according to its intended use.

10.2 Chemical stability

Stable under recommended storage and handling conditions (See section 7).

10.3 Possibility of hazardous reactions

None, when used as directed.

10.4 Conditions to avoid

Heat, naked flames and other ignition sources.

10.5 Incompatible materials

No data available.

Hazardous decomposition products

None if stored, handled and transported properly. In case of fire: see section 5.

SECTION 11: Toxicological information

1

Acute oral toxicity		
No data available		
Acute dermal toxicity		
No data available		
Acute inhalational toxicity		
No data available		
Skin corrosion/irritation		
No Product Name		
1 Steribim® plus (REF 549)	23)	
Source	supplier	
Evaluation/classification	Based on the available data, the classification criteria are not met.	
Serious eye damage/irritation		
No Product Name		
Steribim® plus (REF 549)	23)	
Source	supplier	
Evaluation/classification	Based on available data, the classification criteria are not met.	
Respiratory or skin sensitisat	ion	
No Product Name		
1 Steribim® plus (REF 549)	23)	
Route of exposure	Skin	
Source	supplier	
Evaluation/classification	Based on available data, the classification criteria are not met.	
Route of exposure	respiratory tract	
Source	supplier	
Evaluation/classification	Based on available data, the classification criteria are not met.	
Germ cell mutagenicity		
No Product Name		
1 Steribim® plus (REF 549)	23)	
Source	supplier	
Evaluation/classification	Based on available data, the classification criteria are not met.	
Reproduction toxicity		
No Product Name		
1 Steribim® plus (REF 549)	23)	
Source	supplier	
Evaluation/classification	Based on available data, the classification criteria are not met.	



Trade name: Steribim® plus (REF 54923)

Product no.: 54923

Current version: 1.0.0, issued: 17.01.2022 Replaced version: -, issued: - Region: GB

Card	Carcinogenicity		
No	Product Name		
1	Steribim® plus (REF 54923)		
Source		supplier	
Evaluation/classification		Based on available data, the classification criteria are not met.	

STO	STOT - single exposure		
No	Product Name		
1	Steribim® plus (REF 54923)		
Source		supplier	
Evaluation/classification		Based on available data, the classification criteria are not met.	

STC	STOT - repeated exposure		
No	o Product Name		
1	Steribim® plus (REF 54923)		
Source		supplier	
Evaluation/classification		Based on available data, the classification criteria are not met.	

Aspiration hazard	
Based on available data, the classification criteria are not met.	

11.2 Information on other hazards

Endocrine disrupting properties

No data available.

Other information

No data available.

SECTION 12: Ecological information

12.1 Toxicity

loxicity	to fish ((acute)

No data available

Toxicity to fish (chronic)

No data available

Toxicity to Daphnia (acute)

No data available

Toxicity to Daphnia (chronic)

No data available

Toxicity to algae (acute)

No data available

Toxicity to algae (chronic)

No data available

Bacteria toxicity

No data available

12.2 Persistence and degradability

No data available.

12.3 Bioaccumulative potential

No data available.

12.4 Mobility in soil

No data available.

12.5 Results of PBT and vPvB assessment

Results of PBT and vPvB assessment	
PBT assessment	The components of this product are not considered to be a PBT.



Trade name: Steribim® plus (REF 54923)

Product no.: 54923

Current version: 1.0.0, issued: 17.01.2022 Replaced version: -, issued: - Region: GB

vPvB assessment

The components of this product are not considered to be a vPvB.

12.6 Endocrine disrupting properties

No data available.

12.7 Other adverse effects

No data available.

12.8 Other information

Other information

Do not discharge product unmonitored into the environment.

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Product

Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company. dispose of in accordance with local regulation.

Packaging

Residues must be removed from packaging and when emptied completely disposed of in accordance with the regulations for waste removal. Incompletely emptied packaging must be disposed of in the form of disposal specified by the regional disposer.

SECTION 14: Transport information

14.1 Transport ADR/RID/ADN

The product is not subject to ADR/RID/ADN regulations.

14.2 Transport IMDG

The product is not subject to IMDG regulations.

14.3 Transport ICAO-TI / IATA

The product is not subject to ICAO-TI / IATA regulations.

14.4 Other information

No data available.

14.5 Environmental hazards

Information on environmental hazards, if relevant, please see 14.1 - 14.3.

14.6 Special precautions for user

No data available.

14.7 Maritime transport in bulk according to IMO instruments

Not relevant

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture <u>EU regulations</u>

Regulation (EC) No 1907/2006 (REACH) Annex XIV (List of substances subject to authorisation)

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances considered as substances requiring authorisation as listed on Annex XIV of the REACH regulation (EC) 1907/2006.

REACH candidate list of substances of very high concern (SVHC) for authorisation

According to available data and the information provided by preliminary suppliers, the product does not contain substances that are considered substances meeting the criteria for inclusion in annex XIV (List of Substances Subject to Authorisation) as laid down in Article 57 and article 59 of REACH (EC) 1907/2006.

Regulation (EC) No 1907/2006 (REACH) Annex XVII: RESTRICTIONS ON THE MANUFACTURE, PLACING ON



Trade name: Steribim® plus (REF 54923)

Product no.: 54923

Current version: 1.0.0, issued: 17.01.2022 Replaced version: -, issued: - Region: GB

THE MARKET AND USE OF CERTAIN DANGEROUS SUBSTANCES, MIXTURES AND ARTICLES

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances subject to restriction as listed in Annex XVII of the REACH regulation (EC) 1907/2006.

Directive 2012/18/EU on the control of major-accident hazards involving dangerous substances

This product is not subject to Part 1 or 2 of Annex I.

Other regulations

Adhere to the national sanitary and occupational safety regulations when using this product.

15.2 Chemical safety assessment

A chemical safety assessment has not been carried out for this mixture.

SECTION 16: Other information

Sources of key data used to compile the data sheet:

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

Directives 2000/39/EC, 2006/15/EC, 2009/161/EU, (EU) 2017/164.

National Threshold Limit Values of the corresponding countries as amended in each case.

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.

The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding section.

Creation of the safety data sheet

UMCO GmbH - D-21107 Hamburg, Georg-Wilhelm-Strasse 187, Tel.: +49(40)555 546 300, Fax: +49(40)555 546 357, e-mail: umco@umco.de

This information is based on our present knowledge and experience.

The safety data sheet describes products with a view to safety requirements.

It does not however, constitute a guarantee for any specific product properties and shall not establish a legally valid contractual relationship.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO GmbH.

Prod-ID 782420