

Endo Access Burs

Unique carbide endodontic instruments

MEISINGER'S Endo Access Burs are special carbide endodontic burs for endodontic treatments. The extra-long and very narrow necks, which have a diameter of only 1 mm, allow an outstanding view of the working field and prevent damage from deep root canal walls. Despite this narrow necks the Endo Access Burs have a high rigidity and concentricity. The design of the Endo-drills convince with its high cutting capacity and durable carbide working parts. These burs are heatless and do not create a fine powdery dust.



At a glance

- Special instruments for endodontic treatments
- Long narrow shafts provide dentists with excellent visibility, allowing a better view of deep target areas
- Optimum, vibration-free concentricity
- Shaft design imparts unique effectiveness
- Less dust production than ultrasonic tips
- Color code allows quick identification and selection

Endo Access Burs

	Figure	Shank	Ref.-No.	Size	Total length	Packaged à	Amount	
	191A	RAL (205)	500 205 004 006 005	005	31 mm	2		
	191A	RAL (205)	500 205 004 006 006	006	31 mm	2		
	191A	RAL (205)	500 205 004 006 008	008	31 mm	2		
	191A	RAL (205)	500 205 004 006 010	010	31 mm	2		
	191A	RAL (205)	500 205 004 006 012	012	31 mm	2		
	191A	RAL (205)	500 205 004 006 014	014	31 mm	2		
	191A	RAL (205)	500 205 004 006 018	018	31 mm	2		
	191A	RAXL (206)	500 206 004 006 005	005	34 mm	2		
	191A	RAXL (206)	500 206 004 006 006	006	34 mm	2		
	191A	RAXL (206)	500 206 004 006 008	008	34 mm	2		
	191A	RAXL (206)	500 206 004 006 010	010	34 mm	2		
	191A	RAXL (206)	500 206 004 006 012	012	34 mm	2		
	191A	RAXL (206)	500 206 004 006 014	014	34 mm	2		

CE 0044

Available from the dental trade.

Your trade depot:

Your delivery address:

Hager & Meisinger GmbH

Hansemannstr. 10 • 41468 Neuss • Germany
 Phone: +49 2131 2012-0
 info@meisinger.de • www.meisinger.de

Meisinger USA, L.L.C.

10200 E. Easter Avenue • Centennial, Colorado 80112 • USA
 Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
 Toll free: +1 (866) 634-7464
 info@meisingerusa.com • www.meisingerusa.com

85FL110 - 0915