

Camera drape

- Reduces the risk of infection
- Protects the equipment
- Telescopic folding, easy and quick to drape
- Several different models and sizes
- Flexible models, fit all types of cameras and cables





Telescopic folding

Camera drape

Ref. No.	Dimension (WxL)	Adjustable feature	Insert feature	Dispenser/ Transp. box
705740	8x122 cm	Straight open tip with separate plastic tape and 1 elastic band	None	30/120
705750	15x249 cm	Straight open tip with attached plastic tape and 2 elastic bands	Insert/instruction label	40/80
705820	12x244 cm	Triangle shaped perforated tip with attached plastic tape Integrated closing tape at the end opening	Insert/instruction label	15/60
80004130	14x254 cm	Triangle shaped perforated tip with integrated plastic tape Integrated closing tape at the end opening	Insert/instruction label	15/60
705830	18x246 cm	Triangle shaped elastic tip with attached plastic tape Integrated closing tape at the end opening	Insert/instruction label	15/60
84541211	14x254 cm	Triangle shaped perforated tip with integrated plastic tape	Adjustable cardboard	25/50
705837	16x254 cm	Triangle shaped open tip with attached plastic tape Closing tape attached to the plastic ring	Hard plastic ring with removable cardboard	24/96





Perforated tip

Elastic tip



