








# EMS Handpiece Chart

Ref. No.	Name	Handpiece color	Features	Compatibility	Connector	Light	Thermal Disinfection Ready	Thermal Disinfection validated	DAC validated
EN-034	Original Piezon universal light handpiece	Grey	Light Fiber optic	Light BiK		Yes	No	No	No
EN-037	Original Piezon handpiece PM600	Blue	Metal nozzle	PM600		No	No	No	No
EN-041	Original Piezon universal handpiece	Grey	Universal type 3 pins Black O-ring Plastic nozzle	BiK, miniPiezon, miniMaster, PM400, Kermit, AFS2		No	Yes	No	Yes
EN-046	Original Piezon universal handpiece	Grey	Universal type 3 pins Black O-ring Metal nozzle	BiK, miniPiezon, miniMaster, PM400, Kermit, AFS2		No	Yes	No	Yes
EN-047	Original Piezon handpiece i.Piezon	White	Plastic nozzle	i.Piezon BiK		No	No	No	No
EN-060	Original Piezon LED handpiece	White	New type 4 pins Red O-ring Metal nozzle	BiK NO PAIN, miniMaster LED, PM700, AFMP		Yes	Yes	No	Yes
EN-061	Original Piezon handpiece	White	New type 4 pins Blue O-ring Metal nozzle	BiK NO PAIN, miniMaster LED, PM700, AFMP		No	Yes	No	Yes

