

# VITA PHYSIODENS®

Available in 28 x VITA SYSTEM 3D-MASTER® shades (0M1, 0M3, 1M1 – 5M3) and 15 x VITA classical A1–D4® shades (without B1)

## Upper Anteriors



ovoid

|                                         |                                        |                                       |                                       |
|-----------------------------------------|----------------------------------------|---------------------------------------|---------------------------------------|
| <b>01S*</b><br><br>42.3<br>7.6<br>10.8  | <b>02L</b><br><br>45.7<br>8.4<br>11.5  | <b>03M</b><br><br>47.2<br>8.6<br>10.8 | <b>04L</b><br><br>49.2<br>9.0<br>12.7 |
| <b>05L</b><br><br>50.1<br>9.1<br>11.4   | <b>06L</b><br><br>51.3<br>9.7<br>13.1  |                                       |                                       |
| <b>T6S**</b><br><br>41.0<br>7.4<br>10.4 | <b>T2S*</b><br><br>44.9<br>8.3<br>10.1 | <b>T1S</b><br><br>45.3<br>8.7<br>10.2 | <b>T3M</b><br><br>45.6<br>8.4<br>11.7 |
| <b>T4M*</b><br><br>45.6<br>8.8<br>11.3  | <b>T5M</b><br><br>46.2<br>8.8<br>11.1  | <b>T6L</b><br><br>47.3<br>8.9<br>12.4 | <b>T7M</b><br><br>47.5<br>8.8<br>11.6 |
| <b>T8L*</b><br><br>49.0<br>9.2<br>11.6  | <b>T9L</b><br><br>49.6<br>9.4<br>13.1  |                                       |                                       |



triangular



rectangular

X1M\*

46.3



8.8



11.0

X2L

46.2



8.7



11.6

X3L

47.8



9.0



11.8



square

Z1S

44.1



8.1



9.8

Z2S

47.5



8.7



10.1

## Lower Anteriors

L1S\*

33.0



4.9



8.8

L2M\*

34.5



5.2



10.4

L3M\*

36.3



5.2



10.5

L4M\*

35.5



5.4



9.4

L5M\*

37.0



5.5



10.4

L6L

38.2



5.8



11.4

L7L

38.2



5.5



L8L

43.5



6.3



13.0

\*) Bleached shade(s) 0M1 and 0M3 available exclusively for selected mould(s).

\*\*) Mould(s) only available in VITA SYSTEM 3D-MASTER shade(s).

VITAPAN  
EXCELL

VITAPAN  
PLUS

VITAPAN  
LINGOFORM

VITA  
PHYSIODENS

VITAPAN

VITAPAN  
CUSPIFORM

VITAPAN  
SYNOFORM

VITA  
MFT

VITA VIONIC  
VIGO

VITA VIONIC  
FRAME

|                                 |                                                                                                     |             |             |             |             |             |             |             |             |             |             |             |             |             |             |             |             |
|---------------------------------|-----------------------------------------------------------------------------------------------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| VITA<br>PHYSIODENS<br>Anterior  | UJ / No.<br><br>mm | 01S<br>42.3 | 02L<br>45.7 | 03M<br>47.2 | 04L<br>49.2 | 05L<br>50.1 | 06L<br>51.3 | T6S<br>41.0 | T2S<br>44.9 | T1S<br>45.3 | T3M<br>45.6 | T4M<br>45.6 | T5M<br>46.2 | T6L<br>47.3 | T7M<br>47.5 | T8L<br>49.0 | T9L<br>49.6 |
|                                 | LJ / No.<br><br>mm | L1S<br>33.0 | L2M<br>34.5 | L3M<br>36.3 | L6L<br>38.2 | L6L<br>38.2 | L8L<br>43.5 | L1S<br>33.0 | L2M<br>34.5 | L1S<br>33.0 | L3M<br>36.3 | L4M<br>35.5 | L3M<br>36.3 | L5M<br>37.0 | L7L<br>38.2 | L5M<br>37.0 | L6L<br>38.2 |
| VITA<br>PHYSIODENS<br>Posterior | UJ / No.<br><br>mm | 21E<br>30.5 | 22E<br>32.9 | 23E<br>33.5 | 24E<br>34.8 | 24E<br>34.8 | 25E<br>36.7 | 21E<br>30.5 | 21E<br>30.5 | 20E<br>29.4 | 22E<br>32.9 | 22E<br>32.9 | 23E<br>33.5 | 23E<br>33.5 | 23E<br>33.5 | 24E<br>34.8 | 25E<br>36.7 |
|                                 | LJ / No.<br><br>mm | 21E<br>32.7 | 22E<br>33.7 | 23E<br>34.6 | 24E<br>36.1 | 24E<br>36.1 | 25E<br>37.3 | 21E<br>32.7 | 21E<br>32.7 | 20E<br>31.0 | 22E<br>33.7 | 22E<br>33.7 | 23E<br>34.6 | 23E<br>34.6 | 23E<br>34.6 | 24E<br>36.1 | 25E<br>37.3 |







|                                 |                                                                                                     |             |             |             |             |             |
|---------------------------------|-----------------------------------------------------------------------------------------------------|-------------|-------------|-------------|-------------|-------------|
| VITA<br>PHYSIODENS<br>Anterior  | UJ / No.<br><br>mm | X1M<br>46.3 | X2L<br>46.2 | X3L<br>47.8 | Z1S<br>44.1 | Z2S<br>47.5 |
|                                 | LJ / No.<br><br>mm | L3M<br>36.3 | L4M<br>35.5 | L5M<br>37.0 | L1S<br>33.0 | L5M<br>37.0 |
| VITA<br>PHYSIODENS<br>Posterior | UJ / No.<br><br>mm | 23E<br>33.5 | 23E<br>33.5 | 24E<br>34.8 | 21E<br>30.5 | 22E<br>32.9 |
|                                 | LJ / No.<br><br>mm | 23E<br>34.6 | 23E<br>34.6 | 24E<br>36.1 | 21E<br>32.7 | 22E<br>33.7 |

The information contained in the combination table is only intended as a guide, and can vary depending on the patient's individual situation.

# VITA PHYSIODENS®

Available in 27 x VITA SYSTEM 3D-MASTER® shades (0M3, 1M1 – 5M3) and 15 x VITA classical A1 – D4® shades (without B1)

## Posterior

|                                                                                                                                       |                                                                                                                                        |                                                                                                                                        |
|---------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>20E</b></p> <p>29.4 8.9</p>  <p>31.0 8.3</p>   | <p><b>21E*</b></p> <p>30.5 9.5</p>  <p>32.7 8.7</p>  | <p><b>22E*</b></p> <p>32.9 9.5</p>  <p>33.7 9.0</p> |
| <p><b>23E*</b></p> <p>33.5 11.1</p>  <p>34.6 9.8</p> | <p><b>24E*</b></p> <p>34.8 11.8</p>  <p>36.1 9.8</p> | <p><b>25E</b></p> <p>36.7 11.2</p>  <p>37.3 9.5</p> |

\*) Bleached shade 0M3 available exclusively for selected mould(s).