


Specimen	Length (mm)	Width (mm)	Height (mm)
013	41.7	7.9	9.5
025	42.8	8.1	9.9
034	43.3	8.0	9.7
035	43.2	8.0	9.8
043	43.6	8.2	9.4
086	47.1	8.9	9.9
097	49.2	9.0	11.5
098*	48.1	9.0	11.8
099	52.5	9.8	12.5
T36*	43.7	8.1	10.1
T53	44.6	8.1	9.2
T56	44.4	8.3	10.0
T66	44.6	8.3	10.2
T67	45.1	8.4	10.0
T76*	45.9	8.2	10.4
T77*	46.2	8.5	10.6
T88*	47.1	8.7	10.6



rectangular


square


T98




T99




X13




X66




X77




X87




X96




X99




Z51




Z61




Z74




Z84



Z85



Z97



VITAPAN
EXCELL

VITAPAN
PLUS

VITAPAN
LINGOFORM

VITA
PHYSIODENS

VITAPAN

VITAPAN
CUSPIFORM

VITAPAN
SYNOFORM

VITA
MFT

VITA VIONIC
VIGO

VITA VIONIC
FRAME

Lower Anteriors

Fully anatomical incisal edges	L3 	L5* 	L5L 	L7*
	L9 	L11* 	L13 	L15
Abraded incisal edges	L4* 	L8 	L10 	L12
	L14 			

*) Bleached shade(s) 0M1 and 0M3 available exclusively for selected mould(s).

VITAPAN Anterior	UJ / No.  mm	013 41.7	025 42.8	034 43.3	035 43.2	043 43.6	086 47.1	097 49.2	098 48.1	099 52.5	T36 43.7	T53 44.6	T56 44.4	T66 44.6	T67 45.1	T76 45.9	T77 46.2	T88 47.1	T98 50.5	T99 48.0
	LJ / No.  mm	L4 32.7	L3 33.4	L3 33.4	L3 33.4	L5 35.3	L11 36.9	L8 38.9	L13 37.0	L14 41.7	L4 32.7	L5 35.3	L5 35.3	L5 35.3	L4 32.7	L7 36.2	L11 36.9	L11 36.9	L15 38.0	L11 36.9
VITAPAN CUSPIFORM Posterior	UJ / No.  mm	40C 28.1	40C 28.1	41C 29.9	41C 29.9	41C 29.9	42C 31.7	44C 32.8	42C 31.7	45C 34.8	40C 28.1	41C 29.9	41C 29.9	42C 31.7	42C 31.7	42C 31.7	42C 31.7	44C 32.8	45C 34.8	44C 32.8
	LJ / No.  mm	40C 29.5	40C 29.5	41C 31.3	41C 31.3	41C 31.3	42C 33.4	44C 34.6	42C 33.4	45C 36.3	40C 29.5	41C 31.3	41C 31.3	42C 33.4	42C 33.4	42C 33.4	42C 33.4	44C 34.6	45C 36.3	44C 34.6
VITAPAN SYNOFORM Posterior	UJ / No.  mm	10S 29.0	10S 29.0	11S 30.0	11S 30.0	11S 30.0	12S 31.3	14S 33.1	12S 31.3	14S 33.1	10S 29.0	11S 30.0	11S 30.0	12S 31.3	12S 31.3	12S 31.3	12S 31.3	14S 33.1	14S 33.1	14S 33.1
	LJ / No.  mm	10S 30.5	10S 30.5	11S 31.9	11S 31.9	11S 31.9	12S 32.7	14S 35.4	12S 32.7	14S 35.4	10S 30.5	11S 31.9	11S 31.9	12S 32.7	12S 32.7	12S 32.7	12S 32.7	14S 35.4	14S 35.4	14S 35.4

VITAPAN Anterior	UJ / No.  mm	X13 41.5	X66 45.9	X77 47.0	X87 46.9	X96 50.0	X99 50.6	Z51 44.6	Z61 45.5	Z74 46.7	Z84 46.7	Z85 47.5	Z97 52.2
	LJ / No.  mm	L4 32.7	L9 36.1	L11 36.9	L11 36.9	L10 39.0	L15 38.0	L3 33.4	L3 33.4	L9 36.1	L9 36.1	L13 37.0	L12 41.0
VITAPAN CUSPIFORM Posterior	UJ / No.  mm	41C 29.9	42C 31.7	42C 31.7	44C 32.8	44C 32.8	45C 34.8	43C 31.0	43C 31.0	42C 31.7	42C 31.7	44C 32.8	45C 34.8
	LJ / No.  mm	41C 31.3	42C 33.4	42C 33.4	44C 34.6	44C 34.6	45C 36.3	43C 32.3	43C 32.3	42C 33.4	42C 33.4	44C 34.6	45C 36.3
VITAPAN SYNOFORM Posterior	UJ / No.  mm	11S 30.0	12S 31.3	12S 31.3	14S 33.1	14S 33.1	14S 33.1	13S 31.3	13S 31.3	12S 31.3	12S 31.3	14S 33.1	14S 33.1
	LJ / No.  mm	11S 31.9	12S 32.7	12S 32.7	14S 35.4	14S 35.4	14S 35.4	13S 33.9	13S 33.9	12S 32.7	12S 32.7	14S 35.4	14S 35.4

The information contained in the combination table is only intended as a guide, and can vary depending on the patient's individual situation.