



HU-FRIEDY SIGNATURE SERIES™ PROCEDURE TUBS

GET ORGANIZED
& BE EFFICIENT

How the best perform



IMS CASSETTE AND TUB SYSTEM

HOW CASSETTES AND TUBS WORK TOGETHER:

- Together, procedure tubs and cassettes create the most efficient and organized system for managing instruments and consumable products.
- Time savings achieved by using IMS (Instrument Management System) cassettes and procedure tubs will allow the staff to spend time on more value added activities in the dental office.
- Procedure tubs can be used as portable operator drawers to store, organize and transport consumable materials to and from the storage/sterilization area and the operator.
- All of the materials needed for a specific procedure can be stored in the tub. A different tub and color can be used for each different procedure type (crown & bridge, composite, etc.).
- The Signature Series procedure tubs are available in colors that match IMS cassette rails, which allow the staff to match the tubs and cassettes by procedure.



AN EXAMPLE OF HOW IT WORKS:

- Blue cassette rails and a blue tub may represent a composite procedure. All of the consumable materials for composite procedures can be stored and organized in the blue tub and instrumentation stored in the cassette with blue rails, allowing for quick identification.

PROCEDURE TUB SET-UP EXAMPLES



Hygiene



Sealants



Composite



Endotics



Crown and Bridge

SIGNATURE SERIES PROCEDURE TUB BENEFITS:

STANDARDIZES MATERIAL SET-UP FOR EVERY TYPE OF PROCEDURE

- Increases efficiency of materials management
- Eliminates time-consuming tray preparation for every procedure
- Tubs may be configured to meet your office needs or are available as a complete tub

PROCEDURE TUB COLORS MATCH HU-FRIEDY IMS CASSETTE RAIL COLORS

- Increases office organization and workflow efficiency
- Allows for quick identification of consumable products
- Provides esthetic integration with IMS cassettes

LOCKING COVERS

- Protects items from airborne materials during transportation and storage
- Adheres to infection control standards

MICROBAN® ANTIMICROBIAL PRODUCT PROTECTION

- Microban inhibits the growth of microorganisms that can cause stains, odors and product degradation
- Microban is registered with the EPA for these applications

ALLOWS FOR STORAGE IN CENTRAL STERILIZATION OR EACH OPERATORY

- Tubs can be transported to the operatory prior to each procedure, eliminating inventory management of materials in each operatory

SIGNATURE SERIES™ TUBS

MICROBAN® ANTIMICROBIAL PRODUCT PROTECTION

Microban protection is integrated into the tub and tub components during the manufacturing process. Microban technology works by interrupting key cell functions that allow the microorganisms to function, grow and reproduce. Microban protection begins to work as soon as the microorganism comes into contact with the product surface and works continuously to inhibit the growth of microbes that can cause stain, odors and product degradation - for the lifetime of the product.

Includes:

- Signature Series Tub
- Tub Liner*
- Tub Divider (package of 8)
- Single Tub Cup w/Cover*
- Double Tub Cup w/Cover*
- Divided Slide Tray
- Locking Tub Cover*

SIGNATURE SERIES COMPLETE TUB

12-3/4" x 11-3/16" x 4"

Red	IMS-1400
Gray	IMS-1401
White	IMS-1402
Orange	IMS-1403
Yellow	IMS-1405
Purple	IMS-1406
Blue	IMS-1408
Green	IMS-1409



*Clear components and Tub Liner do not contain Microban®

SIGNATURE SERIES™ TUB ACCESSORIES



SIGNATURE SERIES TUB

12 1/4" x 10 7/8" x 2 3/4"

Red	IMS-1410
Gray	IMS-1411
White	IMS-1412
Orange	IMS-1413
Yellow	IMS-1415
Purple	IMS-1416
Blue	IMS-1418
Green	IMS-1419



LOCKING TUB COVER*

12 3/4" x 11 3/16" x 2 3/4" IMS-1420

- Protects items from airborne materials during transportation and storage



TUB INSERT

9 3/4" x 10 3/4" x 1 1/32" IMS-1423

- The insert features 10 compartments for easy organization of tub materials
- Sits on rim of tub allowing for storage underneath



DIVIDED SLIDE TRAY*

9 3/4" x 5 3/8" x 1 1/16" IMS-1421


- Can be used to store cotton rolls, cotton squares, clamps, syringe tips, etc.



FLAT SLIDE TRAY*

9 3/4" x 5 3/8" x 1 1/16" IMS-1422

- Can be used to store discs, strips, toffemeil, articulating paper, etc.



DIVIDERS

12 3/4" x 11 3/16" x 1 3/8" IMS-1427

- Includes 8 dividers

TUB LINER*

9 5/8" x 8 7/8" x 1/8" IMS-1428

- Non-slip liner prevents small items from slipping under divider sections

SIGNATURE SERIES™ CLIP ON HOLDERS



Capsule Clip-On Holder



Syringe Clip-On Holder

CAPSULE CLIP-ON HOLDER

8" x 2 1/2" x 3/4"

IMS-1429

- Allows quick access to composite material during procedure
- Conveniently clips on any side of tub
- Holds up to 8 capsules, 4 syringe tubes and capsule/gun dispenser

SYRINGE CLIP-ON HOLDER

8" x 2 5/8" x 3/4"

IMS-1430

- Allows quick access to composite material during procedure
- Conveniently clips on any side of tub
- Holds up to 15 syringe tubes and 3 applicator tips

SIGNATURE SERIES™ TUB CUPS AND COVERS

Can be used to store brushes, clamps, cotton pellets, syringe tips, prophy paste, prophy angle tips, pumice, implant tools, etc.



SINGLE TUB WITH COVER*

2 1/4" x 1 5/8" x 1 1/8"

IMS-1424

DOUBLE TUB WITH COVER*

3 3/8" x 2 1/4" x 1 1/8"

IMS-1425

LONG TUB WITH COVER*

3 3/8" x 2 1/4" x 1 1/8"

IMS-1426

- Provides 5 storage compartments
- Dividers are removable for customization
- Perfect for storing capsules, discs, abrasive strips, wedges, matrix bands, etc.

*Clear components and Tub Liner do not contain Microban®

SIGNATURE SERIES™ SYRINGE STANDS



20-UNIT SYRINGE STAND

9" x 4 3/4" x 2 1/2"

IMS-1431

- Stands up for easy access to syringes
- Lays flat to fit inside a tub or drawer
- Holds up to 20 syringe tubes

10-UNIT SYRINGE ADD ON

8 1/2" x 2 3/8" x 3/4"

IMS-1432

- An add-on for the Syringe 20 Unit Stand
- Provides additional space for syringe storage
- Holds up to 10 syringe tubes

SIGNATURE SERIES™ RACKS



Adjustable

ADJUSTABLE 3 TUB RACK

MIN 7 1/2" - MAX 15 1/2" x 10 1/2" x 13"

IMS-1433

- Holds a maximum of 3 tubs with lids or 6 cassettes
- Non-slip feet
- Sturdy metal construction
- Includes fasteners to secure and stabilize
- No screws, tools, nuts or bolts needed

ADJUSTABLE 4 TUB RACK

MIN 7 1/2" - MAX 15 1/2" x 10 1/2" x 18 3/4"

IMS-1434

- Holds a maximum of 4 tubs with lids or 8 cassettes
- Non-slip feet
- Sturdy metal construction
- Includes fasteners to secure and stabilize
- No screws, tools, nuts or bolts needed


ORDERING INFORMATION

Product Description	Part Codes
Complete Tub - Red	IMS-1400
Complete Tub - Gray	IMS-1401
Complete Tub - White	IMS-1402
Complete Tub - Orange	IMS-1403
Complete Tub - Yellow	IMS-1405
Complete Tub - Purple	IMS-1406
Complete Tub - Blue	IMS-1408
Complete Tub - Green	IMS-1409
Tub Only - Red	IMS-1410
Tub Only - Gray	IMS-1411
Tub Only - White	IMS-1412
Tub Only - Orange	IMS-1413
Tub Only - Yellow	IMS-1415
Tub Only - Purple	IMS-1416
Tub Only - Blue	IMS-1418
Tub Only - Green	IMS-1419

Product Description	Part Codes
Locking Tub Cover	IMS-1420
Divided SlideTray	IMS-1421
Flat Slide Tray	IMS-1422
Tub Insert	IMS-1423
Tub Cup Single With Cover	IMS-1424
Tub Cup Double With Cover	IMS-1425
Tub Cup Long With Cover	IMS-1426
Tub Dividers	IMS-1427
Tub Liner	IMS-1428
Capsule Clip-On Holder	IMS-1429
Syringe Clip-On Holder	IMS-1430
20-Unit Syringe Stand	IMS-1431
10-Unit Syringe Add On	IMS-1432
Adjustable 3 Tub Rack	IMS-1433
Adjustable 4 Tub Rack	IMS-1434

(Surface disinfect only. Do not heat sterilize Signature Series Tubs and Accessories.)

For more information about
Hu-Friedy's Total Solutions visit
HU-FRIEDY.COM



Hu-Friedy Mfg. Co., LLC: 3232 N. Rockwell St. | Chicago, IL 60618 | USA
Telephone: 1-800-Hu-Friedy (1-773-975-6100) | **Website:** Hu-Friedy.com

©2012 Hu-Friedy Mfg. Co., LLC. All rights reserved. IMS-380/1012

Hu-Friedy