

let's keep
connected and
feel safe

 **polydentia**
swiss manufacture 

Polydentia SA
via Cantonale 47
6805 Mezzovico-Vira
Switzerland

phone: +41 91 946 29 48
fax: +41 91 946 32 03

info@polydentia.ch
polydentia.ch



Vista-Tec *ilGrande*

The extended protection

Vista-Tec ilGrande is much more than an elegant face protection equipment. It is a statement of your love towards yourself, your practice and your patients. Its wide, resistant and replaceable polycarbonate shield grants a high level of protection against liquid splashes, while the premium optical properties provide pleasant and relaxed vision. Vista-Tec ilGrande also features a headband which makes it comfortable to wear on longer working days. Both frame and headband can be autoclaved for maximum safety.

Vista-Tec ilGrande kit
REF 5663
1 frame, 2 ilGrande shields, 1 headband

Vista-Tec ilGrande shield
REF 5666
5 polycarbonate shields



- extralarge shields for increased protection
- premium optical quality
- CE EN 166:2001 certified
- ultra-light design
- disinfectable and replaceable shields
- autoclavable frame and headband
- optimized comfort and stability



more info at
polydentia.ch

 **polydentia**
swiss manufacture 



Vista-Tec laStoria

The original face shield

The iconic and elegant Vista-Tec face shield, appreciated by dental and health professionals since its conception in 1989, now redesigned with enhanced protection features.

The ultra-lightweight face shield includes a resistant, replaceable polycarbonate shield that provides premium optical properties for a relaxed vision, reducing unnecessary eyestrain over longer work hours.

Vista-Tec laStoria features a silicone headband that enables optimized stability and comfort.

Vista-Tec laStoria kit

REF 5621

1 frame, 2 laStoria shields, 1 headband

Vista-Tec laStoria shield

REF 5626

5 polycarbonate shields

- premium optical quality
- disinfectable and replaceable shields
- ultra-light design
- CE EN 166:2001 certified
- headband for greater comfort and stability



more info at
polydentia.ch

 **polydentia**
swiss manufacture 



Vista-Tec ClipOn

Protective shields for thin-frame spectacles and loupes

ClipOn enables fixing a shield on spectacles. The practical ClipOn clips provide support for both ClipOn shield and the larger L-ClipOn shield, ideal for Galilean loupes. ClipOn clips can be inserted on spectacles or loupes with temple arms of a maximum height of 1cm.

Vista-Tec ClipOn kit

REF 5622

1 pair ClipOn clips, 5 ClipOn shields

Vista-Tec ClipOn shield

REF 5646

10 shields

Vista-Tec L-ClipOn kit

REF 5630

1 pair ClipOn clips, 5 L-ClipOn shields

Vista-Tec L-ClipOn shield

REF 5632

10 shields

