


Intra- and extra-oral.

Appliances for every treatment technique.



DENTAURUM
**ONLINE
SHOP**
shop.dentaurum.com





Complete delivery program for all areas of treatment.

The range of intra-oral and extra-oral products from Dentaureum extends from simple elastics or separators for intra-oral use to innovative appliances such as the FaceFormer.

During the development, production and selection of all intra- and extra-oral products, great care is taken to comply with the Dentaureum quality standards. Not only does this guarantee a high degree of security for the practitioner and patient, but it also enables time-saving and effective treatment.




Elastic auxiliaries

Available in many colors, diameters and with various traction forces.

Dentalastics® separators

For quick and easy separation of teeth.



Extra-oral appliances



Noninium® facebow, medium

nickel-free

(Nickel trace elements of $\leq 0.2\%$ are due to the process of manufacturing)

Break resistant laser welded connection between outer and inner bow.

Safety Note:

Facebows must always be used with headgear and neck bands with safety modules!

Inner bow \emptyset	Safety caps	REF	Quantity
1.15 mm / 45	white	740-000-00	1 piece




Standard facebow, mini

Break resistant laser welded connection between outer and inner bow.

Safety Note:

Facebows must always be used with headgear and neck bands with safety modules!

Inner bow \emptyset	Stop loops	Safety caps	REF	Quantity
1.15 mm / 45		white	740-019-00	1 piece
1.15 mm / 45	97 mm	red	741-049-00	1 piece



Standard facebow, medium

Break resistant laser welded connection between outer and inner bow.

Safety Note:

Facebows must always be used with headgear and neck bands with safety modules!

Inner bow \emptyset	Safety caps	REF	Quantity
1.00 mm / 40	white	740-010-00	1 piece
1.15 mm / 45	white	740-020-00	1 piece


Standard facebow, maxi

Break resistant laser welded connection between outer and inner bow.

Safety Note:

Facebows must always be used with headgear and neck bands with safety modules!

Inner bow ø	Safety caps	REF	Quantity
1.15 mm / 45	white	740-021-00	1 piece



Facebow, medium, with stop loops

Facebow: without coating

Inner bow ø: 1.15 mm / 45

Break resistant laser welded connection between outer and inner bow.

Safety Note:

Facebows must always be used with headgear and neck bands with safety modules!

Outer bow size	Inner bow size	Stop loops	Safety caps / Size	Cuspid hooks	REF	Quantity
49 mm	29 mm	83 mm	yellow / Size 1		741-040-00	1 piece
49 mm	32 mm	90 mm	grey / Size 2		741-045-00	1 piece
49 mm	35 mm	97 mm	red / Size 3		741-050-00	1 piece
49 mm	38 mm	104 mm	green / Size 4		741-055-00	1 piece
49 mm	41 mm	111 mm	blue / Size 5		741-060-00	1 piece
49 mm	35 mm	97 mm	red / Size 3	•	741-150-00	1 piece



Accessories / Extra-oral



Dentalastics® extra-oral plastic chain

Color	Length	REF	Quantity
transparent	118 mm	774-300-00	10 pieces



Extra-oral latex elastics

Highest elasticity and fatigue resistance.

The elastic packs are marked according to size and strength with symbols.

Every type can therefore be quickly identified.

Strength:

Strength ratings are based on the elastic being stretched 3 times the listed inner diameter.

Content:

30 x 100 rings = 3000 pieces or 30 x 50 rings = 1500 pieces



Traction force	Inner ø	REF	Quantity
light 2.2 N (225 g)	6.4 mm / 1/4 inch	772-203-02	3000 pieces
light 2.2 N (225 g)	8.0 mm / 5/16 inch	772-204-02	3000 pieces
light 2.2 N (225 g)	9.5 mm / 3/8 inch	772-205-02	3000 pieces
strong 4.0 N (400 g)	8.0 mm / 5/16 inch	772-208-02	1500 pieces
strong 4.0 N (400 g)	12.7 mm / 1/2 inch	772-210-02	1500 pieces



Lip bumper

With plastic shield and hook for intermaxillary elastics.

Indications:

- for distalization of molars
- for stabilization of molars in the late mixed dentition
- with labial tipping of the anterior teeth caused by tongue pressure

for tubes	REF	Quantity
ø 0.90 mm / 36	748-009-00	1 piece
ø 1.15 mm / 45	748-114-00	1 piece

Stop tube

For facebow and lip bumper.



Inner ø	Length	REF	Quantity
1.15 mm / 45	2 mm	480-914-00	100 pieces

Stop screw

acc. to Ladanyi

For facebow and lip bumper.



Inner ø	REF	Quantity
1.15 mm / 45	480-918-02	10 pieces

Key for hexagon socket screws

For stop screws acc. to Ladanyi.

Only intended for use in the laboratory!



Key width	REF	Quantity
0.9 mm	480-924-02	1 piece

Activator tubes

E.g. for Teuscher activator. These tubes can be built directly into an activator to hold a facebow.

Used for attaching a facebow.



Inner ø	REF	Quantity
1.15 mm / 45	607-300-00	10 pieces
1.30 mm / 51	607-301-00	10 pieces

Elastics and auxiliaries / Intra-oral



Elastic placer

For seating and removing intra-oral elastics.
Color: assorted colors

REF 772-400-00

50 pieces



Intra-oral latex elastics


Highest elasticity and fatigue resistance.
The elastic packs are marked according to size and strength with symbols.
Every type can therefore be quickly identified.

Strength:

Strength ratings are based on the elastic being stretched 3 times the listed inner diameter.

Content:

10 x 100 pieces = 1000 pieces



	Color	Traction force	Inner ø	REF	Quantity
1	natural	light 0.7 N (71 g)	3.2 mm / 1/8 inch	772-103-00	1000 pieces
2	natural	light 0.7 N (71 g)	4.8 mm / 3/16 inch	772-104-00	1000 pieces
3	natural	light 0.7 N (71 g)	6.4 mm / 1/4 inch	772-106-00	1000 pieces
4	natural	light 0.7 N (71 g)	8.0 mm / 5/16 inch	772-107-00	1000 pieces
5	natural	light 0.7 N (71 g)	9.5 mm / 3/8 inch	772-109-00	1000 pieces
6	natural	medium 1.3 N (128 g)	3.2 mm / 1/8 inch	772-113-00	1000 pieces
7	natural	medium 1.3 N (128 g)	4.8 mm / 3/16 inch	772-114-00	1000 pieces
8	natural	medium 1.3 N (128 g)	6.4 mm / 1/4 inch	772-116-00	1000 pieces
9	natural	medium 1.3 N (128 g)	8.0 mm / 5/16 inch	772-117-00	1000 pieces
10	natural	medium 1.3 N (128 g)	9.5 mm / 3/8 inch	772-119-00	1000 pieces
11	natural	strong 1.8 N (184 g)	3.2 mm / 1/8 inch	772-123-00	1000 pieces
12	natural	strong 1.8 N (184 g)	4.8 mm / 3/16 inch	772-124-00	1000 pieces
13	natural	strong 1.8 N (184 g)	6.4 mm / 1/4 inch	772-125-00	1000 pieces
14	natural	strong 1.8 N (184 g)	8.0 mm / 5/16 inch	772-126-00	1000 pieces
15	natural	strong 1.8 N (184 g)	9.5 mm / 3/8 inch	772-127-00	1000 pieces

Intra-oral non-latex elastics

Highest elasticity and fatigue resistance.

The elastic packs are marked according to size and strength with symbols.

Every type can therefore be quickly identified.

Strength:

Strength ratings are based on the elastic being stretched 3 times the listed inner diameter.

Content:

30 x 100 pieces = 3000 pieces



Color	Traction force	Inner ø	REF	Quantity	
transparent	medium 1.3 N (128 g)	3.2 mm / 1/8 inch	772-313-01	3000 pieces	
transparent	medium 1.3 N (128 g)	4.8 mm / 3/16 inch	772-314-01	3000 pieces	
transparent	medium 1.3 N (128 g)	6.4 mm / 1/4 inch	772-316-01	3000 pieces	
transparent	medium 1.3 N (128 g)	8.0 mm / 5/16 inch	772-317-01	3000 pieces	
transparent	strong 1.8 N (184 g)	3.2 mm / 1/8 inch	772-323-01	3000 pieces	
transparent	strong 1.8 N (184 g)	4.8 mm / 3/16 inch	772-324-01	3000 pieces	
transparent	strong 1.8 N (184 g)	6.4 mm / 1/4 inch	772-325-01	3000 pieces	
transparent	strong 1.8 N (184 g)	8.0 mm / 5/16 inch	772-326-01	3000 pieces	



Container for intra-oral elastics

For mini patient packs. With wall mount.

Dimensions (W x H x D): 508 mm x 169 mm x 107 mm

Color: white

REF 772-000-00

1 piece



Intra-oral latex elastics

Special sizes

Cuspid to cuspid, for facebows with cuspid hooks. Highest elasticity and fatigue resistance.

Strength ratings are based on the elastic being stretched 3 times the listed inner diameter.

Traction force	Inner ø	REF	Quantity
medium 1.3 N (128 g)	16 mm / 5/8 inch	772-100-00	1000 pieces
medium 1.3 N (128 g)	19 mm / 3/4 inch	772-200-00	1000 pieces



Elastics and auxiliaries / Intra-oral



Dentalastics® plastic ligatures

For ligating arch wires to brackets, instead of wire ligatures.

Color	Inner ø	Height	REF	Quantity
transparent	1.3 mm	0.9 mm	774-002-01	1000 pieces



Dentalastics® separators

For quick and easy separation of teeth. Radiopaque (ideal setting 60 kV and 3 mA).

Color	Inner ø	Indication	REF	Quantity
green	1.5 mm	for anteriors	774-100-01	1000 pieces
blue	2.1 mm	for posteriors	774-200-01	1000 pieces



SUPER Elasto-Force plastic chains

Consistent traction force over a long period of time. High elasticity with excellent elastic recovery.

Color	Inner ø	Length	Connector	REF	Quantity
transparent	1.3 mm	4.6 m	•	774-220-01	1 roll
pink	1.3 mm	4.6 m	•	774-231-01	1 roll
blue	1.3 mm	4.6 m	•	774-232-01	1 roll
red	1.3 mm	4.6 m	•	774-233-01	1 roll
silver	1.3 mm	4.6 m	•	774-235-01	1 roll
transparent	1.3 mm	4.6 m		774-221-01	1 roll
blue	1.3 mm	4.6 m		774-226-01	1 roll
silver	1.3 mm	4.6 m		774-229-01	1 roll
pearl	1.3 mm	4.6 m		774-236-01	1 roll



Dentalastics® Personal / assortment


Colored ligatures

Especially delicate plastic ligatures, instead of wire ligatures, for ligating archwires in brackets.

Content:

10 colors (without yellow or orange), 5 sticks each, total 50 sticks x 20 rings = 1000 pieces

Color	Inner ø	Height	REF	Quantity
10 colors	1.3 mm	0.7 mm	774-561-04	1 assortment



Dentalastics® Personal / Single packs

Colored ligatures

Especially delicate plastic ligatures, instead of wire ligatures, for ligating archwires in brackets.

Contents: 50 sticks x 20 rings = 1000 pieces

Color	Inner ø	Height	REF	Quantity
silver	1.3 mm	0.7 mm	774-549-01	1 pack
transparent	1.3 mm	0.7 mm	774-551-02	1 pack
black	1.3 mm	0.7 mm	774-552-01	1 pack
white	1.3 mm	0.7 mm	774-553-01	1 pack
pink	1.3 mm	0.7 mm	774-554-01	1 pack
blue	1.3 mm	0.7 mm	774-555-01	1 pack
red	1.3 mm	0.7 mm	774-556-01	1 pack
purple	1.3 mm	0.7 mm	774-557-01	1 pack
yellow	1.3 mm	0.7 mm	774-558-01	1 pack
green	1.3 mm	0.7 mm	774-559-01	1 pack
orange	1.3 mm	0.7 mm	774-560-01	1 pack
pearl	1.3 mm	0.7 mm	774-563-01	1 pack



Oral screens



Oral screen Ulmer model

Indications:

- incorrect functions of the tongue (pressing of the tongue, incorrect position of the tongue)
- biting the lips
- sucking the lips
- thumb sucking
- weak muscle structure of the lips

Advantages:

- no pressure points due to the flexible material
- more comfortable to wear
- better acceptance by the patient
- easy to monitor the performance of the myofunctional exercises

Special appliance containers see Supplies!

Size	REF	Quantity
Mini / Size 1 = for preschool age patients	074-000-00	1 piece
Maxi / Size 2 = for patients from school age to adults	074-001-00	1 piece



FaceFormer ONE blue with hygiene box

acc. to Berndsen

FaceFormer therapy targets the cause and is holistic. Regular FaceFormer training actively tones tissues, muscles and skin in the mouth and face. It strengthens and harmonizes muscle functions that are particularly important for chewing, swallowing, speaking, breathing, tooth position as well as facial expressions and tightening of the skin.

Indications:

- difficulty in swallowing (dysphagia)
- incorrect functions of the tongue
- dysfunctional mastication
- speech disorders with functional causes
- curing thumb/pacifier sucking, lip biting etc.
- paralysis of the oral and facial areas
- following an operation on the face, jaw or tongue
- malalignment of the teeth, anomalies in the jaw
- grinding the teeth (bruxism), temporomandibular problems
- sleep disorders, snoring, sleep apnea
- hearing loss, tinnitus
- dizziness, Menière's disease
- problems compensating pressure in the middle ear
- activation of nasal breathing, curing mouth breathing
- etc.

Color	Form	REF	Quantity
transparent blue	anatomical	074-101-01	1 piece

