

TRANSCEND™

UNIVERSAL COMPOSITE

Restorations with Just One Shade

Transcend composite is designed with Ultradent's Resin Particle Match™ technology which allows the Universal Body shade to blend with almost any tooth color anywhere in the mouth, no blocker needed. If you prefer a layering technique, or for more complex anterior cases, Transcend composite also includes four dentin shades and two enamel shades.

Before

After



Deep staining from amalgam tattooing presents one of the most difficult restoration situations to clinicians. In this case only the Transcend composite Universal Body shade was used to replace the amalgam, no blocker needed. Note the excellent color blending of the preserved oblique ridge.



Transcend composite allows you to complete restorations with just one shade. Thanks to Ultradent's Resin Particle Match™ technology, the refraction indices of the resin and particles work together to allow Transcend composite to blend with the surrounding tooth color. That means you can use Transcend composite Universal Body shade almost anywhere in the mouth and know that it will look natural and beautiful, even in larger restorations. Plus, Transcend composite features ideal handling for manipulation and sculptability.



Transcend Syringe 1pk
4 g syringe

Dentin	1pk	Enamel	1pk
A1D	4727	EN	4731
A2D	4728	EW	4732
A3D	4729	Body	1pk
B1D	4730	UB	4733



Transcend Singles 1pk
10 x 0.2 g singles

Dentin	1pk	Enamel	1pk
A1D	4744	EN	4748
A2D	4745	EW	4749
A3D	4746	Body	1pk
B1D	4747	UB	4757



4734 - Transcend UB Syringe 4pk
4 x 4 g syringes Universal Body shade



4817 - Transcend UB Singles 4pk
10 x 0.2 g singles Universal Body shade



4726 - Transcend Syringe Intro Kit
1 x 4 g syringe of each shade:
A1D, A2D, A3D, B1D, EN, EW, UB



4814 - Transcend Singles Intro Kit
10 x 0.2 g singles of each shade:
A1D, A2D, A3D, B1D, EN, EW, UB

Pair Composite Wetting Resin with any Ultradent composite to improve instrument and composite glide when sculpting and contouring.



3059 - Composite Wetting Resin Syringe 2pk
2 x 1.2 ml syringes

